







Click domes and click sleeves selection guide.



Life sounds brilliant.

Products

Click Dome Types

Products	Description	Size
	The open dome is for patients who require maximum comfort and no occlusion.	4 mm 8 mm 10 mm
	The semi-open dome combines the benefits of an open and a closed dome for patients who desire the comfort of an open fitting but need more power.	8-12 mm
	The closed dome is for patients who require more loudness than the open dome can offer.	6 mm 8 mm 10 mm
	The double dome is for patients who require more loudness and more retention in the ear.	8 mm 10 mm
	The closed click sleeve is for patients who require more comfort and greater retention.	XS S M L
	The vented click sleeve combines the benefit of greater retention and the desire for the comfort of an open fitting.	XS S M L

Measurement instructions:

Place the measuring tool vertically, against the ear. To determine the correct receiver size, find the center of the ear canal and read directly across to the corresponding number (see arrow). If the reading is in between numbers, err on the side of the larger number and longer receiver length.



PROPER placement



IMPROPER placement

RIC length measuring tool.

Cut out the measuring tools included here and hook them over your ear for an accurate measurement.

