

Programmierverbindungen Januar 2021



		Active Active Pro	Styletto X	Pure Charge&Go X	CROS Pure Charge&Go X	Pure 312 X	CROS Pure 312 X	Motion Charge&Go X	Motion Charge&Go PX	Motion Charge&Go SPX	Silk X	Cros Silk X	Styletto Connet	Styletto 2/1	Styletto	Pure 10 Nx	Pure 312 Nx	CROS Pure 312 Nx	Pure 13 Nx	Pure Charge&Go Nx	CROS Pure Charge&Go Nx	Motion 13 Nx	Motion 13P Nx	Motion Charge&Go Nx	Silk Nx	CROS Silk Nx	Insio Nx IIC	Insio Nx CIC/ITC/ITE ohne BT	Insio Nx ITC/ITE mit BT	Intuis 3 RIC 312	Intuis 3 S	Intuis 3 M	Intuis 3 P/SP	Prompt S	Prompt P/SP			
KABELGEBUNDENE PROGRAMMIERUNG	Adapter Typ 10	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
	Adapter Typ 312	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	•	-	-	-	-	-	-	-	-	-	-	-	•	•	-	-	•	-		
	Adapter Typ 13	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	-	-	-	•	•	-	-	-	-	-	-	-	-	•	-	-	-	-	
	Flex Connector	-	-	-	-	-	-	-	-	-	•	•	-	-	-	-	-	-	-	-	-	-	-	-	•	•	-	•	•	-	-	-	-	-	-	-	-	-
	Programmierbuchse CS44	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	Programmier-schuh + Kabel	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	-	•	-
WIRELESS-PROGRAMMIERUNG	Noahlink Wireless	•	•	•	•	•	•	•	•	•	-	-	•	-	-	-	•	•	•	•	•	•	•	•	-	-	-	-	•	-	-	-	-	-	-	-	-	
	ConnexAir	-	-	-	-	-	-	-	-	-	•	•	-	•	•	-	-	-	-	-	-	-	-	-	•	•	•	•*	-	-	-	-	-	-	-	-	-	-

*Nur in Verbindung mit e2e wireless Technologie.

Die aufgeführten Informationen können ohne vorherige Ankündigung geändert werden.
(Zuletzt aktualisiert: Januar 2021)