

About Signia

Signia is one of the leading hearing aid brands in the world. Our goal is to enhance human performance through iconic innovations that shape the prescription audio category. We refuse to accept hearing loss as a limitation. It's not about correcting a loss, it's about gaining an edge. Not just performing, performing brilliantly. Not just hearing again, but hearing better than ever. Hearing every word, beat, and breath. With Signia, you can Be Brilliant.

Find out more on signia-pro.com

Or ask your hearing care professional for more information on Signia hearing aids, useful accessories and the Signia app.



Compatible with Android and iOS. The app is free of charge.

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by WSAUD A/S is under license. Other trademarks and trade names are those of their respective owners.

The 'Qi' symbol is a trademark of the Wireless Power Consortium.

** Data on file

*** Study conducted at the University of Northern Colorado, 2015, examined the effectiveness of new Signia features by collecting and analyzing ongoing EEG data while subjects performed speech testing. For both Signia features SpeechMaster and EchoShield, the objective brain behavior measures revealed a significant reduction in listening effort when the feature was activated.

Legal Manufacturer
WSAUD A/S
Nymøllevej 6
3540 Lynge
Denmark

Order No. 04401-99C1-7600 | Printed in Germany
©10.2020 WSAUD A/S

Find us on

 signia-pro.com

 facebook.com/signiahearing

 youtube.com/signiahearing

 instagram.com/signiahearing

Be Brilliant™

Product Portfolio 2021



Be Brilliant™

Signia is one of the leading hearing aid brands in the world. Our goal is to enhance human performance through iconic innovations that shape the prescription audio category. We refuse to accept hearing loss as a limitation. It's not about correcting a loss, it's about gaining an edge. Not just performing, performing brilliantly. Not just hearing again, but hearing better than ever. Hearing every word, beat, and breath.

With Signia, you can Be Brilliant.

Contents




Xperience platform	4
Active Pro	6
Motion X	7
Styletto X	10
Silk X	11
Pure X	12
Styletto Connect / Styletto	14
Nx platform	16
Essential segment	26
Basic segment	32
App and accessories	38
Hearing care professional support	42
Color overview	44
Detailed feature overview	46

Signia Xperience

Empowering people to hear what matters

Performance levels



How do you convince people with hearing loss who have not yet decided to wear hearing aids to do so? At Signia, we have dedicated ourselves to this challenge by finding ways for our newest solutions to help their wearers Be Brilliant. People with hearing loss do not want to feel old because they aren't old. They want hearing solutions that match their active lifestyle and can-do spirit.

Our latest hearing innovations offer leading audiological performance in iconic form factors. They combine our proven Signia Xperience YourSound technology built on revolutionary acoustic-motion sensors with Li-ion rechargeability and Bluetooth connectivity so your clients always hear what matters, with ease.



Xperience

Signia Xperience is available across all performance levels from 7 to 1 in four additional devices, so you can offer all your customers outstanding hearing performance.

Click ITE

BTE



Active Pro



Motion C&G X



Motion C&G PX



Motion C&G SP X

SLIM-RIC

Click ITE

RIC



Styletto X



Silk X



Pure Charge&Go X



Pure 312 X

Signia Xperience feature overview

Our revolutionary platform Signia Xperience gives hearing aids a new sense to truly understand what the wearer needs to hear even in noise, all combined with advanced connectivity and leading remote assistance.

Features	Information	Performance levels						
		7X	5X	3X	2X	1X	Active Pro	Active
YourSound Technology	Redesigned acoustic sensors plus an integrated motion sensor provide detailed information about where the wearers are and what they are doing. Dynamic Soundscape Processing and Own Voice Processing (OVP™) process this information delivering the most clear and personalized sound.	●●●●●	●●●●●	●●●●●	●●●○	●●●○	●●●●●	●●●○
Acoustic Sensors	Detects more variables of sound than ever before, providing more accurate information about the environment.	✓	✓	✓	✓	✓	✓	✓
Dynamic Soundscape Processing	Steers Sound Clarity and Sound Quality features for natural sound and speech in every situation from any direction, even when moving.	●●●●●	●●●●●	●●●●●	●●●○	●●●○	●●●●●	●●●○
Motion Sensor^A	Detects whether the wearer is walking, running or not. This information is used for Dynamic Soundscape Processing.	✓	✓	✓				
OVP™ (Own Voice Processing)^A	Utilizes real-time recognition of the wearer's voice to deliver a natural own-voice impression.	●●●●●	●●●●●	●●●●●				
Signal Processing / Gain & MPO	Channels/handles in Connexx.	48/20 ^B	32/16	24/12	16/8	16/8	48/20	16/8
Signia Assistant^A	Offers the wearers instant 24/7 support and troubleshooting around the clock. It provides detailed information in Connexx to address each client's fine-tuning needs more precisely.	✓	✓	✓	✓	✓	✓	✓
TeleCare	Provides virtual appointments, remote fine-tuning, progress updates, hearing lessons and daily satisfaction ratings.	✓	✓	✓	✓	✓	✓	✓
Tinnitus	A choice of tinnitus treatments based on amplified notch therapy or traditional noise therapy signals. It is even possible to combine the two approaches.	✓	✓	✓	✓			✓
eWindScreen™	Reduces the annoyance of wind noise when outdoors, improving sound comfort. eWindScreen Binaural requires an e2e wireless link.	✓ Binaural Auto in Universal program	✓ Binaural manual	✓ Monaural	✓ Monaural			✓ Binaural Auto in Universal program
Binaural OneMic Directionality 2.0	Delivers a highly sophisticated speech-in-noise performance. Available for Silk X.	✓	✓	✓				✓
Spatial SpeechFocus	Detects the loudest voice in the environment and adjusts the directional beam towards that voice (i.e. front, side, or backward facing). Available for bilateral fittings with directional microphones and e2e wireless link.	✓ front / back / right / left ^C	✓ front / back / omni ^D					✓ front / back / right / left ^C
Extended Bandwidth	12 kHz bandwidth for enhanced processing of high-frequency speech and environmental sounds.	✓						✓
EchoShield***	Dedicated program for reverberant environments.	✓						✓

●●●●● = highest feature performance
 ●●●● = mid-range feature performance
 ●●● = base feature performance

A) Not available with Silk X
 B) For Motion C&G SP X: 34/16
 C) In Universal program: automatic mic beam to the right / left / rear directions
 D) In Universal program: automatic mic beam to the rear, surround mode for side speech. In manual programs, microphone beam to the right / left / rear directions

Prescription audio that performs

New

DYNAMIC SOUNDSCAPE PROCESSING
DSP

ULTRA HD
E2E



Key features

- Iconic earbud design bridging the gap to high tech consumer products
- Signia Xperience technology for brilliant sound and speech understanding
- Remote control
- TeleCare

Signia app ready

Instant Fit

Equipped with Bluetooth

Recharge-on-the-go



Charging case
Signia Active Pro comes in a pocket-sized Qi wireless charging case that lets wearers put the buds in and out of their ears when they need them.

Active Pro

Performance levels

Active Pro

Active

Our game-changing Signia Active Pro takes hearing aids to that next level for you: Prescription audio that enhances human performance and enables wearers to Be Brilliant. Signia Active can be put in and taken out of the ears as easily as any earbud. Modern-looking, easy to use and fast to fit, these handy devices mark a whole new category of hearing solutions – for everyone who doesn't want to wear hearing aids but wants to enjoy hearing everything that's important to them.

ITE

Model	
Performance	
Matrix (2 ccm)	42/118 dB
Battery	
Type	Lithium-ion power cell

Experience the full potential of BTEs

New

ACOUSTIC MOTION SENSORS
AMS

DYNAMIC SOUNDSCAPE PROCESSING
DSP

OWN VOICE PROCESSING
OVP

ULTRA HD
E2E



Key features

- Most natural own voice and best speech understanding in noise, even when the wearer is in motion
- World's first integrated acoustic-motion sensors
- Long-lasting high quality audio streaming
- TeleCare
- Remote Control
- Connectivity & Streaming with StreamLine Mic and StreamLine TV

Signia app ready

Smaller Design

Equipped with Bluetooth

Rechargeable



Charging case
The Inductive Charger II comes with a charger lid to protect the hearing aids, a fast charging option, and a drying function. Its smart design accommodates custom molds.

Motion Charge&Go X

Performance levels

7X

5X

3X

2X

1X

sDemo

Comparable in its small size to a 312- device, the new Motion Charge&Go X is the ultimate BTE combining sleek design with rechargeability, Bluetooth, and the Signia Assistant. It is flexible to fit, easy to use and comfortable to wear, making this an easy choice for those looking for convenience and robustness – all in one.

BTE

Model		
Performance		
Matrix (2 ccm)	60 dB/133 dB SPL	60 dB/125 dB SPL
Battery		
Type	Lithium-ion power cell	Lithium-ion power cell

Hear what really matters

New

ACOUSTIC MOTION SENSORS **AMS**

DYNAMIC SOUNDSCAPE PROCESSING **DSP**

OWN VOICE PROCESSING **OVP**

ULTRA HD **E2E**



Performance levels

- 7X
- 5X
- 3X
- 2X
- 1X
- sDemo

Motion Charge&Go P X

Smaller than its predecessor without compromising on power, the new Motion Charge&Go P X combines all the latest Signia Xperience hearing technology with the complete convenience of Li-ion recharging, making annoying battery changes a thing of the past. This powerful hearing solution for every situation with telecoil ensures that whatever you spend your day doing, you will be able to perform at your best.

Key features

- Most natural own voice and best speech understanding in noise, even when the wearer is in motion
- World's first integrated acoustic-motion sensors
- Long-lasting high quality audio streaming
- TeleCare
- Remote Control
- Connectivity & Streaming with StreamLine Mic and StreamLine TV

- Signia app ready
- Smaller Design
- Equipped with Bluetooth
- Rechargeable
- Integrated T-Coil



Charging cases
The consumer can select between two charger options, a standard charger or a Dry&Clean charger.

BTE

Model	Earhook	ThinTube 3.0
Performance	77 dB/135 dB SPL	66 dB/130 dB SPL
Battery	Lithium-ion power cell	
Type	Lithium-ion power cell	

Superpower your hearing

New

ACOUSTIC MOTION SENSORS **AMS**

DYNAMIC SOUNDSCAPE PROCESSING **DSP**

OWN VOICE PROCESSING **OVP**

ULTRA HD **E2E**



Performance levels

- 7X
- 5X
- 3X
- 2X
- 1X
- sDemo

Motion Charge&Go SP X

This is the world's first rechargeable Super Power BTE. Closing the gap between Super Power and super features, the innovative Motion Charge&Go SP X ultra-power device is the complete package for those who want to Be Brilliant despite severe to profound hearing loss. With Bluetooth connectivity, Li-ion rechargeability and telecoil, it offers a world of new hearing and handling possibilities not familiar to those with conventional SP devices.

Key features

- Most natural own voice and best speech understanding in noise, even when the wearer is in motion
- World's first integrated acoustic-motion sensors
- Long-lasting high quality audio streaming
- TeleCare
- Remote Control
- Connectivity & Streaming with StreamLine Mic and StreamLine TV

- Signia app ready
- Smaller Design
- Equipped with Bluetooth
- Rechargeable
- Integrated T-Coil




Charging cases
The consumer can select between two charger options, a standard charger or a Dry&Clean charger.

BTE

Model	Motion C&G SP X
Performance	82 dB/140 dB SPL
Battery	Lithium-ion power cell
Type	Lithium-ion power cell

Hearing without compromise



Key features

- Uniquely stylish, award-winning design*
- Signia Xperience technology with the world's first acoustic-motion sensors for clear speech comprehension even in noisy environments
- Tinnitus options
- Remote Control
- Connectivity and high quality audio streaming

- Signia app ready
- Rechargeable
- Equipped with Bluetooth



Charging case
The world's first hearing aid charger with Qi technology – wearers simply place the charger on a charging pad instead of having to plug it into the mains.

Performance levels

- 7X
- 5X
- 3X
- 2X
- 1X
- sDemo

Styletto X

Styletto X is the complete package hearing aid ideal for anyone who doesn't want to compromise on their active lifestyle just because of hearing loss. The world's first SLIM-RIC combines brilliant, award-winning design* and portable rechargeability with exchangeable receivers, the clear sound and speech understanding of our revolutionary Signia Xperience platform, and Bluetooth connectivity. This unique form factor can significantly increase your conversion of people with hearing loss into satisfied clients. When consumers were given a selection of hearing aids to choose from, Styletto increased conversion by 18%.


SLIM-RIC

Model	SLIM-RIC		
Receiver type	S-Receiver	M-Receiver	P-Receiver
Performance			
Matrix (2 ccm)	45/108 dB SPL	60/119 dB SPL	70 /124 dB SPL
Battery			
Type	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell

Source: Hakvoort, Burton: Increasing Style, Reducing Stigma: The Styletto Solution (Signia Backgrounder, 2018)

*Red Dot Design Award 2019 for Styletto Connect, iF Design Award 2019 for Styletto


Discretion guaranteed



Key features

- Nearly invisible in the ear for complete discretion
- Signia Xperience technology for a superior hearing experience
- Click Sleeves for a quick and easy instant-fit
- Remote control
- TeleCare

- Signia app ready
- Compatible with CROS device for people with unaidable hearing loss on one side



One click and you're done.
The sleeves come in four different sizes and a closed and open build version suitable for most ear anatomies and hearing needs.

Performance levels

- 7X
- 5X
- 3X
- 2X
- 1X
- sDemo


Silk X

Silk X is the ideal solution for those who value discretion above all. The only ready-to-wear CIC available on the Signia Xperience platform delivers clear speech understanding in an ultra-discreet size with instant fit. It is tailor-made without the wait: Your clients can walk in, experience Silk X and walk out the door with a new CIC without having to wait for a custom-built device. Thanks to its flexible silicone Click Sleeves, Silk X sits securely in the ear for a brilliant fit on the spot.

Click ITE

Model	Click ITE
Receiver type	Click CIC
Performance	
Matrix (2 ccm)	50/114 dB SPL
Battery	
Type	
Size	10
Power consumption	1.2 mA

Supreme style and streaming with charging-on-the-go



Performance levels

- 7X
- 5X
- 3X
- 2X
- 1X
- sDemo

Pure Charge&Go X

Pure Charge&Go X combines all the advantages of Signia Xperience with Li-ion rechargeability and fully-featured Bluetooth connectivity. Its 20% extra battery capacity provides full flexibility and convenience for higher wearer satisfaction, with streaming of music, calls and TV in high-quality sound. All in a remarkably small design. Simply Charge&Go!

Key features

- Most natural own voice and best speech understanding in noise, even when the wearer is in motion
- World's first integrated acoustic-motion sensors
- Long-lasting high quality audio streaming
- TeleCare
- Remote Control
- Connectivity & Streaming with StreamLine Mic and StreamLine TV

- Signia app ready
- Rechargeable
- Equipped with Bluetooth
- Compatible with CROS device for people with unaidable hearing loss on one side




Charging case
The Inductive Charger II comes with a charger lid to protect the hearing aids, a fast charging option, and a drying function. Its smart design accommodates custom molds.

RIC

Model	RIC			
Receiver type	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
Performance				
Matrix (2 ccm)	45/108 dB	60/119 dB	70/124 dB	75/130 dB
Battery				
Type	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell

Fully versatile, fully connected



Performance levels

- 7X
- 5X
- 3X
- 2X
- 1X
- sDemo

Pure 312 X

Equipped with all the benefits of the Signia Xperience platform, Pure 312 X attracts wearers with its elegantly slim design. Despite its discreet form, Pure 312 X provides fully-featured Bluetooth connectivity for maximum enjoyment of music, calls and TV. An integrated T-coil is also available.

Key features

- Most natural own voice and best speech understanding in noise, even when the wearer is in motion
- World's first integrated acoustic-motion sensors
- Long-lasting high quality audio streaming
- TeleCare
- Remote Control
- Connectivity & Streaming

- Signia app ready
- Equipped with Bluetooth
- Optional T-Coil
- Compatible with CROS device for people with unaidable hearing loss on one side

RIC

Model	RIC			
Receiver type	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
Performance				
Matrix (2 ccm)	45/108 dB	60/119dB	70/124 dB	75/130 dB
Battery				
Size	312	312	312	312
Power consumption	1.7 mA	1.9 mA	1.8 mA	1.8 mA

Supreme style and streaming with charging-on-the-go



Key features

- Stylishly designed new form factor
- Most natural own voice and best speech understanding in noise thanks to Nx chip
- TeleCare
- Remote Control
- Connectivity & Streaming

- Signia app ready
- Rechargeable
- Equipped with Bluetooth



Charging case
Wearers enjoy full-day wearing and charging-on-the-go via Signia Styletto's ultra-slim portable charging case.

Performance levels

- 7
- 5
- 3
- sDemo

Styletto Connect

Signia's new Styletto Connect combines all the best of audiology, convenience and design. Powered by Signia's groundbreaking Nx platform, this fashionable hearwear combines lithium-ion power cells with all benefits of Bluetooth connectivity for recharging and streaming-on-the-go. In addition to its high-class performance, Styletto Connect comes in an exceptional and distinctive new design – with a unique appeal that will attract even more customers.

SLIM-RIC

Model	
Receiver type	miniReceiver 2.0
	
Performance	
Matrix (2 ccm)	60/119 dB
Battery	
Type	Lithium-ion power cell

Embracing life with stylish high-tech



Key features

- Stylishly designed form factor
- Outstanding speech understanding in noise
- TeleCare
- Remote Control

- Signia app ready
- Rechargeable



Charging case
Wearers enjoy full-day wearing and charging-on-the-go via Signia Styletto's ultra-slim portable charging case.

Performance levels

- 7
- 5
- 3
- 2
- 1
- sDemo

Styletto

Transforming the image of hearing aids, Signia's revolutionary hearwear Styletto combines most of the audiological advantages of our acclaimed Signia Nx platform with a stylishly designed form factor. And now it is accessible to even more people: Available in all performance levels, it provides the latest advancements in terms of energy-efficient chip design, unique lithium-ion power cells, and natural sound. All in a refined slim-line styling.

SLIM-RIC

Model	
Receiver type	miniReceiver 2.0
	
Performance	
Matrix (2 ccm)	60/119 dB
Battery	
Type	Lithium-ion power cell

Signia Nx

Outstanding soundscape processing

Performance levels

- 7Nx
- 5Nx
- 3Nx
- 2Nx
- 1Nx

The sophisticated chip technology of Signia Nx packs even more performance into even less space. Its leading technology enables advanced connectivity in the smallest sized hearing aid with Ultra HD e2e and Bluetooth.



Signia Nx

Pure Charge&Go Nx, Pure 312 Nx, 13 Nx and 10 Nx hearing aids together with the Motion Charge&Go Nx, Motion 13 Nx, 13P Nx and Silk Nx cover all performance ranges from performance levels 7 to 1, expanding the range of devices available with Nx technology. That means more choice for your customers, more opportunities for you.

Signia Nx is available across all performance levels from 7 to 1 in six additional devices, so you can offer all your customers outstanding hearing performance.

RIC



Pure Charge&Go Nx



Pure 312 Nx



Pure 13 Nx



Pure 10 Nx

BTE



Motion Charge&Go Nx



Motion 13 Nx



Motion 13P Nx

Click ITE



Silk Nx

ITE



Insio Nx

Signia Nx feature overview

Our game-changing platform Signia Nx delivers the most natural own voice and speech understanding in noise, while combining the most advanced connectivity systems.

Features	Information	Performance levels
		7Nx 5Nx 3Nx 2Nx 1Nx
Audiology		
OVP™ (Own Voice Processing)	Own Voice Processing (OVP™) utilizes real-time recognition of the wearer's voice to deliver a natural own-voice impression. Available for SLIM-RIC/RIC/BTE devices with Ultra HD e2e.	●●●●● ●●●●● ●●●●● - -
3D Classifier	Unlike conventional classification systems which rely on acoustic data, 3D Classifier analyzes multi-dimensional information from the wearer's environment and communication behavior. Information on the wearer's voice activity, motion state and acoustic data are used to tailor the hearing aids' response to the wearer's needs. Motion detection available for Bluetooth-enabled devices.	●●●●● ●●●●● ●●● - -
Sound Clarity	Sound Clarity provides lifelike sound which covers the full spectrum of wearing situations. Exceptional clarity is achieved with features such as: HD Spatial, extended dynamic range, extended bandwidth, EchoShield, HD Music, eWindScreen™, speech and noise management and feedback cancellation.	●●●●● ●●●●● ●●● ●● ●●
Speech Quality	Speech Quality delivers improved understanding of other people's voices in noise by using an exclusive technology called Ultra HD e2e. Features which utilize Ultra HD e2e include narrow directionality, Spatial SpeechFocus™ and TwinPhone.	●●●●● ●●●●● ●●● ●● ●●
Connectivity		
Direct Streaming	Signia Nx hearing aids can be connected directly to Apple devices for phone calls and audio streaming – without additional accessories. With the StreamLine Mic accessory, direct streaming is possible also from Android devices and other Bluetooth mobile devices. With the StreamLine TV accessory, high-quality stereo sound can be streamed directly to Nx hearing aids. Available for Bluetooth-equipped devices.	●●●●● ●●●●● ●●●●● ●●●●● ●●●●●
You Adjust	Technically inclined wearers often appreciate having extra control over their hearing aids. Signia Nx hearing aids offer these wearers a comprehensive range of adjustment options via Signia app.	●●●●● ●●●●● ●●●●● ●●●●● ●●●●●
Remote hearing care		
TeleCare	TeleCare provides remote services including additional tools to follow-up the patients within the whole journey, beginning with pre-appointment and even continuing after purchase. It contributes to a more successful trial phase and ensures higher flexibility. In addition, the user engagement and autonomy during this whole journey will be increased via the Signia app.	●●●●● ●●●●● ●●●●● ●●●●● ●●●●●

●●●●● = highest feature performance

Convenient wireless rechargeability and Bluetooth on all levels



Key features

- Rechargeable lithium-ion power cell
- Long-lasting high quality audio streaming
- TeleCare
- Remote Control
- Connectivity & Streaming

- Signia app ready
- Rechargeable
- Equipped with Bluetooth



StreamLine TV
With StreamLine TV, the TV signal is directly streamed into Pure Charge&Go Nx hearing aids.

StreamLine Mic
The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any device with Bluetooth connectivity.

Performance levels

- 7Nx
- 5Nx
- 3Nx
- 2Nx
- 1Nx
- sDemo

Pure Charge&Go Nx

Simply Charge&Go: Now, Signia's popular Pure Charge&Go Nx with convenient wireless lithium-ion rechargeability, is accessible for an even broader range of users. Also available in all performance levels, it provides top Bluetooth connectivity, remote support with Signia TeleCare and the most high-tech wireless charging in a small design. The proven Nx Fit fitting formula will let wearers feel its simplicity right from the start!

RIC

Model	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
Receiver type				
Performance	45/108 dB	60/119 dB	70/124 dB	75/130 dB
Battery	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell

Small and fully-featured with direct streaming



Key features

- Most natural own voice and best speech understanding in noise thanks to Nx chip
- Small, discreet design
- TeleCare
- Remote Control
- Connectivity & Streaming

- Signia app ready
- Equipped with Bluetooth
- Compatible with CROS device for people with unaidable hearing loss on one side



StreamLine TV
With StreamLine TV, the TV signal is directly streamed into Pure 312 Nx hearing aids.

StreamLine Mic
The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any device with Bluetooth connectivity.

Performance levels

- 7Nx
- 5Nx
- 3Nx
- 2Nx
- 1Nx
- sDemo

Pure 312 Nx

Equipped with all benefits of the Signia Nx platform, Pure 312 Nx keeps the wearer's own voice and remaining soundscape in perfect balance. This way, it offers the most natural own voice and best speech understanding in noise in a small, elegant housing. Despite its discreet design, Pure 312 Nx provides supreme connectivity via Bluetooth-based direct streaming and the Signia app*. Plus, TeleCare lets your customers enjoy full live remote support and fast, smooth trials.

* For Android and iOS.

RIC

Model	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
Receiver type				
Performance	45/108 dB	60/119dB	70/124 dB	75/130 dB
Battery	312	312	312	312
Size	1.2 mA	1.4 mA	1.3 mA	1.3 mA
Power consumption				

Superior connectivity combined with the longest streaming time



Key features

- Most natural own voice and best speech understanding in noise thanks to Nx chip
- Longest wearing time while streaming
- TeleCare
- Remote Control
- Connectivity & Streaming

- Signia app ready
- Equipped with Bluetooth
- Standard battery door can be exchanged for one with an integrated T-coil
- Compatible with CROS device for people with unaidable hearing loss on one side



StreamLine TV
With StreamLine TV, the TV signal is directly streamed into Pure 13 Nx hearing aids.

StreamLine Mic
The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any device with Bluetooth connectivity.

Pure 13 Nx

By separately processing the wearer's own voice and integrating it smoothly into the soundscape, Signia's Pure 13 Nx provides the most natural own voice and best speech understanding in noise. Thanks to Bluetooth, it comes with superior connectivity with direct streaming and the Signia app*. Users enjoy the longest wearing time while streaming. Full live remote support via TeleCare quickly familiarizes them with their Pure 13 Nx hearing aids.

* For Android and iOS.

RIC

Model	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
Receiver type				
Performance				
Matrix (2 ccm)	45/108 dB	60/119 dB	70/124 dB	75/130 dB
Battery				
Size	13	13	13	13
Power consumption	1.2 mA	1.4 mA	1.3 mA	1.3 mA

Simplicity with a naturally clear sound

Key features

- Ultra-small design
- Outstanding speech understanding in noise
- Remote control option
- TeleCare
- Remote Control



Signia app and miniPocket
To enhance its discretion, Pure 10 Nx is remotely controllable via the Signia app and the miniPocket.

Performance levels

- 7Nx
- 5Nx
- 3Nx
- 2Nx
- 1Nx
- sDemo

Pure 10 Nx

The ultra-small Pure 10 Nx expertly combines a highly discreet design with a more natural sound. Its sleek appearance is complemented by a wealth of audiological innovations that produce a clear sound and excellent speech understanding in noise. The non-wireless Pure 10 Nx is the ideal solution for those who wish their hearing aids to be simple and discreet.

RIC

Model	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
Receiver type				
Performance				
Matrix (2 ccm)	45/108 dB	55/113 dB	60/118 dB	65/123 dB
Battery				
Size	10	10	10	10
Power consumption	1.2 mA	1.3 mA	1.2 mA	1.2 mA

The ultimate rechargeable BTE with Bluetooth for almost all needs



Key features

- Rechargeable lithium-ion power cell
- Long-lasting high quality audio streaming
- TeleCare
- Remote Control
- Connectivity & Streaming

- Signia app ready
- Rechargeable
- Equipped with Bluetooth



StreamLine TV
With StreamLine TV, the TV signal is directly streamed into Motion Charge&Go Nx hearing aids.

StreamLine Mic
The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any device with Bluetooth connectivity.

Performance levels

- 7Nx
- 5Nx
- 3Nx
- 2Nx
- 1Nx
- sDemo

Motion Charge&Go Nx

With Motion Charge&Go Nx, Signia's all-round BTE just got even more versatile: Now also available in all performance levels, these models combine all the advantages of Signia Nx with Li-ion powered wireless rechargeability and seamless Bluetooth connectivity in a dependable design. Intuitive to handle, their outstanding sound quality, convenient streaming options, and easy fittings make Motion Charge&Go Nx hearing aids the optimum solution for almost any customer.

BTE

Model	Motion Charge&Go Nx	Motion Charge&Go Nx
Performance		
Matrix (2 ccm)	ThinTube 56/126 dB	Earhook 63/134 dB
Battery		
Type	Lithium-ion power cell	Lithium-ion power cell

Direct streaming with the most versatile fitting options



Key features

- Most natural own voice and best speech understanding in noise thanks to Nx chip
- For individual handling and fitting requirements
- TeleCare
- Remote Control
- Connectivity & Streaming

- Signia app ready
- Equipped with Bluetooth
- Standard battery door can be exchanged for one with an integrated T-coil
- Compatible with CROS device for people with unaidable hearing loss on one side



StreamLine TV
With StreamLine TV, the TV signal is directly streamed into Motion 13 Nx hearing aids.

StreamLine Mic
The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any device with Bluetooth connectivity.

Performance levels

- 7Nx
- 5Nx
- 3Nx
- 2Nx
- 1Nx
- sDemo

Motion 13 Nx

Available in a small model and a sleek P model, the incredibly versatile Motion 13 Nx BTE is ideal for wearers with individual fitting requirements and moderate to severe or even profound hearing loss. It harmonizes the wearer's voice and surrounding sounds, providing the most natural own voice and best speech understanding in noise. Motion 13 Nx and Motion 13P Nx provide full Bluetooth connectivity and Signia TeleCare with an optional modular T-coil integrated in the battery door.

BTE



Model	Motion 13 Nx	Motion 13P Nx
Performance		
Matrix (2 ccm)	ThinTube 56/126 dB Earhook 63/134 dB	ThinTube 65/130 dB Earhook 77/135 dB
Battery		
Size	13	13
Power consumption	1.4 mA	1.7 - 2.0 mA

A practically invisible natural experience. Simply fit and go



Key features

- Best speech understanding in noise and a most natural sound thanks to Nx chip
- Click Sleeves for a quick and easy fit
- Even smaller than its predecessor
- Invisible and discreet
- Binaural OneMic directionality
- TeleCare
- Remote Control



One click and you're done. Available in four sizes, Click Sleeves simply click on to Silk Nx hearing aids and sit securely in the ear.

Performance levels

- 7Nx
- 5Nx
- 3Nx
- 2Nx
- 1Nx
- sDemo

Silk Nx

The Silk™ Nx delivers a natural hearing experience in its discreet size, instant fit and exceptional sound quality. The only ready-to-wear CIC with Ultra HD e2e for superior wireless binaural processing. Thanks to its flexible silicone Click Sleeves, Silk Nx offers a precise fit without waiting times. And for wearers with unaidable hearing loss in one ear, CROS Silk Nx provides the most discreet solution – the world's first CROS CIC hearing aid.

Click ITE



Model	Click CIC
Performance	50/113 dB
Matrix (2 ccm)	50/113 dB
Battery	
Size	10
Power consumption	1.3 mA

The most discreet natural hearing experience with full connectivity



Key features

- Best speech understanding in noise and a most natural sound thanks to Nx chip
- Small size up to high power matrix
- Full wireless and remote control options for all models
- TeleCare
- Remote Control
- Connectivity & Streaming*



StreamLine TV With StreamLine TV, the TV signal is directly streamed into Bluetooth-enabled Inzio Nx ITC and ITE hearing aids.

StreamLine Mic The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any Bluetooth-enabled Inzio Nx ITC/ITEs.

Performance levels

- 7Nx
- 5Nx
- 3Nx
- 2Nx
- 1Nx

Inzio Nx

Signia's Inzio™ Nx ITC and ITE models provide top Bluetooth connectivity with best-in-class sound quality. Tailor-made for the most discerning wearers, Inzio Nx delivers a more natural sound with maximum discretion. Equipped with Ultra HD e2e for superior wireless binaural processing, including as a small CIC with a single microphone, Inzio offers crystal clear sound and unrivalled speech understanding even in very noisy environments.

ITE

Model	IIC	CIC	ITC	ITE
Performance				
Matrix (2 ccm)	50/113 dB1	S 50/113 dB1 M 55/118 dB1 P 65/124 dB1	S 50/113 dB1 M 55/118 dB1 P 65/124 dB1	M 55/118 dB1 P 65/124 dB1
Battery				
Size	10	10	312	312
Power consumption	1.3 mA1	1.3 mA1	1.3 mA1	1.3 mA1

(1) Max. values

* Available with Bluetooth-enabled Inzio Nx ITC and ITE hearing aids

Essential

Proven technology, outstanding comfort

Performance levels

Essential

With the Intuis 3 family from Signia, you can provide your customers with the convenience of easy, discreet handling options and TeleCare. For the first time ever, this generation of Intuis hearing aids – from RIC and BTEs to ITEs to Click ITEs – is remotely controllable via the Signia app and even TeleCare ready.

All RIC, BTE, Click ITE and ITE models in Signia’s Intuis 3 family offer proven technology with 6 channel MPO for advanced fitting options and for more individual adjustments. The P model comes in a new ergonomic design for outstanding comfort. Most models are equipped with automatic directional microphones that provide a more precise hearing experience.

RIC



Intuis 3 RIC312

BTE



Intuis 3 BTEs

ITE



Intuis 3 Click ITEs and ITEs

Essential segment feature overview

Proven technology, outstanding comfort: In a broad range of hearing aid models, the Essential line offers a reliable hearing experience, a high level of wearing comfort and easy handling for almost every customer.

Features

Better hearing

Frequency channels / Compression controls
Number of programs
Directional microphones A
Frequency compression

Sound comfort

Feedback cancellation
Speech and noise management
eWindScreen
e2e wireless 3.0
e2e wireless 2.0
SoundSmoothing

Fitting / Individuality

Data logging
Acclimatization manager
Tinnitus therapy

Audio streamer / Hearing aid control

Signia app
Smart Connexx app and Signia app
Mini Pocket

A) Not for CIC B) Only BTE

Performance levels

Essential

Intuis 3

12/6

4

automatic B

•

•

on/off

-

-

-

-

•

-

-

•

-

•

Small and powerful

Key features

- For mild to severe hearing loss
- Small design for a high level of wearing comfort
- TeleCare
- Remote Control



Signia app and miniPocket
To enhance its discretion, Intuis 3 RIC312 is remotely controllable via the Signia app and the miniPocket.

Intuis 3 RIC312

Performance levels

Essential

Welcome to the family! Completing the proven Signia Intuis 3 range, the remotely controllable Intuis 3 RIC312 is an attractive and convenient solution for mild to severe hearing loss. Thanks to its automatic directionality and a variety of noise reduction features, it offers a reliable hearing experience and improved speech intelligibility even in noisy environments. TeleCare lets your customers enjoy live remote support and fast, smooth trials.

RIC

Model	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
Receiver type				
Performance	45/108 dB	60/119 dB	70/124 dB	75/130 dB
Matrix (2 ccm)				
Battery				
Type	312	312	312	312
Power consumption	0.9 mA	1.0 mA	1.0 mA	1.1 mA

Proven technology, outstanding comfort

Key features

- For all levels of hearing loss
- Ergonomic design for a high level of wearing comfort
- Remotely controllable for a maximum of convenience and discretion
- TeleCare
- Remote Control



Intuis 3 BTEs

Performance levels

Essential

All Signia's Intuis 3 BTE models offer proven technology with 6 channel MPO for the most advanced fitting options and for more individual adjustments. Suitable for all levels of hearing loss, the easy to handle Intuis 3 BTEs are equipped with automatic directional microphones for a precise hearing experience, offer remote control options, and are TeleCare ready for smoother, faster trials and higher customer satisfaction.

BTE

Model	Intuis 3 S	Intuis 3 M	Intuis 3 P	Intuis 3 SP
Receiver type				
Performance	LifeTube 45/124 dB Earhook 55/124 dB	LifeTube 53/125 dB Earhook 60/130 dB	Earhook 70/134 dB	Earhook 80/140 dB
Matrix (2 ccm)				
Battery				
Type	312	13	13	675
Power consumption	0.9 mA	1.0 mA	1.4 mA	2.4 mA

The perfect match

Key features

- For mild to severe hearing loss
- 6 channel MPO for the most advanced fitting options and more individual adjustments
- Remotely controllable for a maximum of discretion and convenience
- TeleCare
- Remote Control



Intuis 3 Click ITEs and ITEs

Performance levels


Essential

With Intuis 3 Click ITEs and tailor-made ITEs, you can offer your customers the full range of in-the-ear models: No matter if they wish for a solution that fits right away or prefer a customized model, Intuis 3 is the perfect match. Equipped with proven technology, all models offer a precise hearing experience, are remotely controllable for a maximum of discretion and TeleCare ready for successful trials.

Click ITE/ITE

Model	Intuis 3 Click CIC, Click ITC	Intuis 3 CIC, ITC, ITE
		
Performance	Click CIC 50/113 dB Click ITC 55/116 dB	CIC 65/124 dB1 ITC 65/124 dB1 ITE 65/124 dB1
Battery		
Size	Click CIC 10, Click ITC 312	CIC 10, ITC 312, ITE 13
Power consumption	1.1 mA	1.1 mA1

(1) Max. values



Basic

Help the world hear better

Basic

The Basic segment of hearing aids offers fitting opportunities that enable you to bring the promise of better hearing to people who would otherwise have no hope of a solution. Easy to learn and easy to fit, even in mobile situations, you can take better hearing to where it's needed – even to outlying areas. This range of BTE and Click ITE devices is suitable for all levels of hearing loss, sizes of ear canal and individual preference.

BTE



Click ITE



ITE



Basic segment feature overview

Available in a wide range of fittings and suitable for all levels of hearing loss from mild to profound, the Basic range offers a good hearing experience for every customer.

Features	Prompt						
	S	P	SP	Click ITEs	CIC	ITC	ITE
Matrix	LifeTube 45/124 dB Earhook 55/124 dB	70/134 dB	80/140 dB	Click CIC 50/113 dB Click ITC 55/116 dB	65/124 dB1	65/124 dB1	65/124 dB1
Channels/Controls	8/4	8/4	8/4	8/4	8/4	8/4	8/4
Programs	4	4	4	Click CIC: 1 Click ITC: 4	4	4	4
Directional microphones	●	●	●	-	-	-	-
T-coil	●	●	●	-	-	○	○
Audio input	-	●	●	-	-	-	-
Noise reduction	situation based	situation based	situation based	situation based	situation based	situation based	situation based
SoundBalance™	●	●	●	●	●	●	●
Feedback cancellation	●	●	●	●	●	●	●

Features	Run				Fun		Fast
	S	P	SP	Click ITEs	P	SP	P
Matrix	LifeTube 45/124 dB Earhook 55/124 dB	70/132 dB	79/140 dB	Click CIC 50/113 dB Click ITC 55/116 dB	70/132 dB	79/140 dB	62/130 dB
Channels	8	8	8	8	6	6	4
Programs	6 preconfigured sound profiles	6 preconfigured sound profiles	6 preconfigured sound profiles	6 preconfigured sound profiles	4 preconfigured sound profiles	4 preconfigured sound profiles	3 preconfigured sound profiles
Directional microphones	●	●	-	-	-	-	-
T-coil	-	-	-	-	-	-	-
Audio input	-	-	-	-	-	-	-
Noise reduction	situation based	situation based	situation based	situation based	always on	always on	always on
SoundBalance™	●	●	●	●	-	-	-
Feedback cancellation	●	●	●	●	●	●	●

(1) Max. values

Make your customers the center of attention

Key features

- A full range of BTE, Click ITE and ITE models
- For all levels of hearing loss
- Remotely controllable for a maximum of convenience and discretion
- TeleCare
- Remote Control



Basic The Prompt family

Signia's Prompt family offers outstanding wearing comfort and ease of use for every form of hearing loss. All models offer remote control options via the Signia app and the miniPocket™ and are a convenient solution for almost all your customers.



BTE

Model	Prompt S	Prompt P	Prompt SP
Performance			
Matrix (2 ccm)	ThinTube 45/124 dB Earhook 55/124 dB	Earhook 70/134 dB	Earhook 80/140 dB
Battery			
Size	312	13	675
Power consumption	0.9 mA	1.4 mA	2.4 mA

Click ITE/ITE

Model	Prompt Click CIC, Click ITC	Prompt CIC, ITC, ITE
Performance		
Matrix (2 ccm)	Click CIC 50/113 dB Click ITC 55/116 dB	CIC 65/124 dB1 ITC 65/124 dB1 ITE 65/124 dB1
Battery		
Size	Click CIC 10, Click ITC 312	CIC 10, ITC 312/10, ITE 13/312
Power consumption	1.1 mA	1.1 mA1

(1) Max. values

Take hearing to a new level

Key features

- Quick and easy fittings everywhere via the MobileFitting app
- For mild to profound hearing loss
- Ergonomic shape for a high level of wearing comfort
- TeleCare
- Remote Control



Basic Run

Easy to fit in just a few minutes with the convenient MobileFitting app, Run from Signia is ideal for on-the-spot fittings. Incredibly comfortable thanks to its ergonomic shape, the convenient and easy-to-use Run delivers an outstanding hearing experience.



BTE

Model	Run S	Run P	Run SP
Performance			
Matrix (2 ccm)	ThinTube 45/124 dB Earhook 55/124 dB	Earhook 70/132 dB	Earhook 79/140 dB
Battery			
Size	312	13	675
Power consumption	0.9 mA	1.3 mA	2.2 mA

Click ITE

Model	Run Click CIC	Run Click ITC
Performance		
Matrix (2 ccm)	50/113 dB	55/116 dB
Battery		
Size	10	312
Power consumption	1.1 mA	1.1 mA

Provide an outstanding hearing experience

Key features

- A high level of amplification for severe to profound hearing loss
- Ergonomic shape for a high level of wearing comfort
- Simple fittings everywhere via the rocker switch

Deliver a powerful listening experience

Key features

- Simple fittings everywhere via the rocker switch
- For mild to severe hearing loss
- Easy handling

Basic

Fun

Equipped with Signia's proven hearing aid technology, Fun covers all essential customer requirements. It offers superb sound quality for severe to profound hearing loss, excellent wearing comfort and is simply fitted by pushing the handy rocker switch.



Basic

Fast

Fast BTE hearing aids offer a high level of performance and deliver just the right amount of amplification for customers with severe hearing loss. They are easily adjusted via the rocker switch for simplified fitting. Sleek and comfortable with a range of features, Signia's Fast ensures an outstanding hearing experience suitable for every customer.



BTE

Model	Fun P	Fun SP
		
Performance		
Matrix (2 ccm)	Earhook 70/132 dB	Earhook 79/140 dB
Battery		
Size	13	675
Power consumption	1.2 mA	2.2 mA

BTE

Model	Fast P
	
Performance	
Matrix (2 ccm)	Earhook 62/130 dB
Battery	
Size	13
Power consumption	1.1 mA

App, accessories and hearing aid compatibility

This overview shows you which apps and accessories work with our hearing aids. Depending on the solution your customers choose, they can even enhance their hearing experience and convenience with this well-matching portfolio.

Hearing aid	Signia app	CROS/ BiCROS	Streamers and transmitters		Remote controls	Charging stations	Inductive Charger II	Standard Charger	Dry&Clean Charger
		CROS	StreamLine TV	StreamLine Mic	miniPocket	Charger	Charger	Charger	Charger
Active Pro	●	-	●	●	●	●	-	-	-
Motion Charge&Go X	●	●	●	●	●	●	●	-	-
Motion Charge&Go P X	●	●	●	●	●	●	-	●	●
Motion Charge&Go SP X	●	●	●	●	●	●	-	●	●
Styletto X	●	-	●	●	●	●	-	-	-
Silk X	●	●	-	-	●	-	-	-	-
Pure Charge&Go X	●	●	●	●	●	●	●	-	-
Pure 312 X	●	●	●	●	●	-	-	-	-
Styletto Connect	●	-	●	●	●	● ^A	-	-	-
Styletto	●	-	-	-	●	● ^B	-	-	-
Pure 10 Nx	●	-	-	-	●	-	-	-	-
Pure 312 Nx	●	● ^C	●	●	●	-	-	-	-
Pure 13 Nx	●	● ^C	●	●	●	-	-	-	-
Pure Charge&Go Nx	●	● ^C	●	●	●	-	● ^D	-	-
Motion 13 Nx	●	● ^C	●	●	●	-	-	-	-
Motion 13P Nx	●	● ^C	●	●	●	-	-	-	-
Motion Charge&Go Nx	●	● ^C	●	●	●	-	● ^D	-	-
Silk Nx	●	● ^E	-	-	●	-	-	-	-
Insio Nx	●	● ^{G,H}	● ^F	● ^F	●	-	-	-	-
Intuis 3	●	-	-	-	●	-	-	-	-
Intuis 3 RIC312	●	-	-	-	●	-	-	-	-
Prompt	●	-	-	-	●	-	-	-	-
Run	●	-	-	-	●	-	-	-	-

A) Compatible with Styletto Connect charging case
 B) Compatible with Styletto charging case
 C) Compatible with CROS Pure 312 Nx or CROS Pure Charge&Go Nx
 D) Compatible with the Pure Charge&Go/Motion Charge&Go/Cellion charger

E) Compatible with CROS Silk Nx
 F) Compatible with Bluetooth-enabled Insio Nx ITC/ITE
 G) Insio Nx IIC and CIC are compatible with CROS Silk Nx
 H) Bluetooth-enabled Insio Nx ITC and ITE are compatible with CROS Pure 312 Nx

Accessory overview

Streamers



Via StreamLine TV, TV sound can be fed directly into Bluetooth-equipped Signia hearing aids. Wearers can even adjust the TV volume individually via the Signia app.



The StreamLine Mic transforms Signia devices with Bluetooth into more than state-of-the-art hearing aids. They also become hands-free mobile headsets and wireless stereo headphones that stream highest quality audio into both ears from any Bluetooth device. Wearers can even adjust the streaming volume individually via the Signia app.

Remote controls



Small enough to carry on a keyring, the miniPocket is the ideal solution for anyone who wants discreet, convenient hearing aid control without the need for a large remote control or smartphone.

CROS/BiCROS



For people with unaidable hearing loss on one side, Signia CROS devices receive sound on that side and transmit it to a hearing aid on the other side.

Charging stations



The chargers recharge hearing aids so your clients don't need to worry about handling batteries. The Inductive Charger II for Motion C&G X, Pure Charge&Go X, Pure Charge&Go Nx and Motion Charge&Go Nx also dries hearing aids while charging to give a full day's ultra-durable hearing performance.



The new Dry & Clean Charger for Motion C&G P X and Motion C&G SP X offers a UV-C cleaning functionality while charging and drying as an additional benefit, while the Standard Charger gives your clients the option of a basic economic charging alternative without a dry and clean function.

All your clients' needs, all together

Signia app Greater convenience

Simpler all round, more convenient for all, the new Signia app combines all Signia's existing apps into one unified environment for all the user's needs. And because it includes TeleCare, it integrates seamlessly into your fitting process, making it easier for your customers to stay in touch with you when they are away from your premises. Suitable for use with all Signia hearing aids, its smart user interface only shows the user functions applicable to their own devices.



Compatible with Android and iOS.
The app is free of charge.

TeleCare

- Stay in closer contact with your customers:
- Real-time contact via text message, voice call or video call
- Intuitive monitoring of your customers' daily satisfaction
- Set virtual appointments for the precise remote adjustment of the hearing aids

Remote Control

- Personalize your customers' hearing experience further:
- Hearing program and volume control
- Sound Balance Adjustment
- 360° Spatial Configurator for individual adjustment of directionality
- Continuous monitoring of the hearing aids' connection and battery status

Connectivity & Streaming

- Easy management of audio streaming accessories (e.g. StreamLine TV, StreamLine Mic) as well as CROS devices
- Volume level control while streaming
- Adaptive streaming volume



TeleCare



Remote Control



Connectivity & Streaming



Digital Assistant



Signia Assistant Your very own hearing companion


Imagine being able to help your clients to Be Brilliant every single step of their hearing journey outside your office? Now you can – with the revolutionary Signia Assistant.

Get closer to your clients

Signia Assistant is the smart new way to better address your clients' needs. It acts as an extension of yourself so you can see the challenges they experience in their daily lives and support them in optimizing their hearing experience. All they need is their Signia Xperience hearing aids and the Signia app.

See every detail with Connexx

With the Connexx fitting software, you'll find all the changes made by the Assistant. It's the smart solution for more satisfied customers.



With the Signia Assistant you can:

- Give your clients the assurance of 24/7 support
- Increase your client satisfaction and conversion rate with the power to address hearing aid wearers' needs more accurately
- Increase the quality of follow up appointments through a new level of real-world usage insights
- Fine-tune with ease by seeing exactly how your clients choose to personalize their hearing experience



Perfection is in the detail
Components for hearing aids



From RIC miniReceivers™ to BTE earpieces and coupling options to Click Sleeves 2.0 for Click ITEs and RICs, we supply you with all necessary items so you can tailor the perfect hearing solution for your customers.

Hearing aid components

RICs*: miniReceivers S, M, P, HP
 Click Sleeves 2.0 XS, S, M, L
 Click Molds™

BTEs: Thin Tubes
 Ear Tips
 Basic Coupling
 Earmolds

Click ITEs: Click Sleeves 2.0 XS, S, M, L

* Styletto is only compatible to M-size receivers and all sizes of Click Sleeves.

Meeting all your needs

The Signia hearing care professional support

Every day you try to meet and exceed your customers' hearing needs and wishes. In Signia you have a partner who fully supports you. Our broad range of hearing aid components, counseling and trial support as well as fitting tools enhances your workflow and provides you with all you need to find the best individual solution for your customers.

For detailed information, see our website signia-pro.com

Reducing effort, increasing satisfaction

Support for counseling, fittings and trials



Signia offers you a broad range of counseling and fitting support – from programming software and tools to mobile fitting units to hearing aids for demonstration and advisory tools – up to comprehensive remote hearing care. This enhances your workflow and makes trials and fittings smoother and more pleasant for your customers.

Tools and support for counseling, fittings and trials

Connex 9
 iScan™ II
 Unity™ 3
 ConnexAir
 Noahlink
 Noahlink Wireless
 HI-PRO 2 USB
 TeleCare
 SoundDemo Toolkit

Detailed feature overview SLIM-RIC and RICs

Type	Styletto X	Pure Charge&Go X	Pure 312 X	Styletto Connect	Styletto
Performance levels	7X/5X/3X/2X/1X			7/5/3	7/5/3/2/1
Picture					
Product details					
Matrix (2 ccm)	S 45/108 dB M 60/119 dB P 70/124 dB -	S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB	S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB	- 60/119 dB -	- 60/119 dB -
Battery size	Lithium-ion power cell	Lithium-ion power cell	312	Lithium-ion power cell	Lithium-ion power cell
Power consumption			1.7 – 1.9 mA		
Features					
Frequency channels/ Compression controls	7X – 48/20 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8	7X – 48/20 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8	7X – 48/20 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8	7 – 48/20 5 – 32/16 3 – 24/12 -	7 – 48/20 5 – 32/16 3 – 24/12 2 – 16/8 1 – 16/8
Hearing programs	6	6	6	6	6
e2e wireless	●	● ¹³	● ¹³	● ¹³	● ¹³
Direct streaming	●	●	●	●	-
TwinPhone	-	● ¹⁶	● ¹⁶	● ¹⁵	● ¹⁵
TeleCare	●	● ¹⁹	● ¹⁹	● ¹⁹	● ¹⁸
Tinnitus therapies (available on 7/5/3/2 PL)	●	●	●	●	●
CROS/BiCROS	-	●	●	-	-
Push button	-	-	-	-	-
Volume control	-	● ²⁴	● ²⁴	-	-
Rechargeable	●	●	-	●	●
T-coil	-	-	○	-	-
Battery door on/off function	-	-	●	-	-
Audio input	● ²⁸	● ²⁸	● ²⁸	● ²⁸	-
Modifications					
Color conversion kit	-	○	○	-	-
App/Accessory compatibility					
Signia app	●	●	●	●	●
Smart Connect	-	-	-	-	-
StreamLine TV/StreamLine Mic	●	●	●	●	-
miniPocket	●	●	●	●	●
Connexx Smart Remote	-	-	-	-	-
ConnexxAir	-	-	-	-	●
Noahlink Wireless	●	●	●	●	-
Programming adapter/cable	-	-	-	-	-
Fitting set	miniReceiver 2.0	miniReceiver 2.0	miniReceiver 2.0	-	-

Pure Charge&Go Nx	Pure 312 Nx	Pure 13 Nx	Pure 10 Nx	Intuis 3 RIC312
7Nx/5Nx/3Nx/2Nx/1Nx				Essential
				
S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB	S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB	S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB	S 45/108 dB M 55/113 dB P 60/118 dB HP 65/123 dB	S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB
Lithium-ion power cell	312	13	10	312
	1.2 – 1.4 mA	1.2 – 1.4 mA	1.2 – 1.3 mA	0.9 – 1.1 mA
7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	12/6
6	6	6	6	4
● ¹³	● ¹³	● ¹³	-	-
●	●	●	-	-
● ¹⁷	● ¹⁶	● ¹⁶	-	-
● ¹⁹	● ¹⁹	● ¹⁹	● ¹⁸	● ¹⁸
●	●	●	●	-
● ^{2,21}	● ^{2,21}	● ^{2,21}	-	-
●	-	-	○	● ²³
● ²⁴	● ²⁴	● ²⁴	-	● ²³
●	-	-	-	-
-	-	○	-	●
-	●	●	●	●
● ²⁸	● ²⁸	● ²⁸	-	-
○	○	○	○	○
●	●	●	●	●
-	-	-	-	-
●	●	●	●	●
-	-	-	-	-
●	●	●	-	-
-	size 312	size 13	size 10	size 312
miniReceiver 2.0				

(2) Not available with PL 1
 (13) Ultra HD e2e with binaural audio signal processing
 (15) To be activated via Signia app or miniPocket
 (16) To be activated via rocker switch, Signia app or miniPocket
 (17) To be activated via push button, Signia app or miniPocket
 (18) Basic remote tuning via 4-band-equalizer

(19) Full live remote tuning via Connexx
 (21) With either CROS Pure 312 Nx or CROS Pure Charge&Go Nx
 (23) Exchangeable controls – programmable via Connexx
 (24) Programmable via Connexx
 (28) Via direct streaming

● Standard
 ○ Optional

Detailed feature overview BTEs

Type	Motion C&G X	Motion C&G P X	Motion C&G SP X	Motion C&G Nx	Motion 13 Nx	
Performance levels	7X/5X/3X/2X/1X			7Nx/5Nx/3Nx/2Nx/1Nx		
Picture						
Product details						
Matrix (2 ccm)	Earhook ThinTube/ LifeTube 3.0	60 dB/133 dB SPL 60 dB/125 dB SPL -	77 dB/135 dB SPL 66 dB/130 dB SPL -	82 dB/140 dB SPL -	56/126 dB - 63/134 dB	56/126 dB - 63/134 dB
Battery size	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	13	
Power consumption					1,4 – 2,0 mA	
Features						
Frequency channels/ Compression controls	7X – 48/20 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8	7X – 48/20 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8	7X – 34/16 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	
Hearing programs	6	6	6	6	6	
e2e wireless	●	●	●	● ¹³	● ¹³	
Direct streaming	●	●	●	●	●	
TwinPhone	-	-	-	● ¹⁷	● ¹⁶	
TeleCare	●	●	●	● ¹⁹	● ¹⁹	
Tinnitus therapies	-	-	-	●	●	
CROS/BiCROS	●	●	●	● ^{2,21}	● ^{2,21}	
Push button	-	-	-	●	-	
Volume control	-	-	-	● ²⁴	● ²⁴	
Rechargeable	●	●	●	●	-	
T-coil	-	●	●	-	○	
Battery door on/off function	-	-	-	-	●	
Audio input	-	-	-	● ²⁸	● ²⁸	
Modifications						
Color conversion kit	○	○	○	○	○	
App/Accessory compatibility						
Signia app	●	●	●	●	●	
Smart Connect	-	-	-	-	-	
StreamLine TV/StreamLine Mic	●	●	●	●	●	
miniPocket	●	●	●	●	●	
easyPocket	-	-	-	-	-	
ConnexxAir	-	-	-	-	-	
Noahlink Wireless	●	●	●	●	●	
Programming adapter/cable	-	-	-	-	size 13	
Fitting set	ThinTube/Earhook	ThinTube/Earhook	Earhook	ThinTube/Earhook	ThinTube/Earhook	

Motion 13P Nx	Intuis 3 S	Intuis 3 M	Intuis 3 P	Intuis 3 SP	
	Essential				
Performance levels	65/130 dB - 77/135 dB	- 45/124 dB 55/124 dB	- 53/125 dB 60/130 dB	- 70/134 dB	- 80/140 dB
Battery size	13	312	13	13	675
Power consumption	1,7 – 2,0 mA	0,9 mA	1,0 mA	1,4 mA	2,4 mA
Features					
Frequency channels/ Compression controls	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	12/6	12/6	12/6	12/6
Hearing programs	6	4	4	4	4
e2e wireless	● ¹³	-	-	-	-
Direct streaming	●	-	-	-	-
TwinPhone	● ¹⁶	-	-	-	-
TeleCare	● ¹⁹	● ¹⁸	● ¹⁸	● ¹⁸	● ¹⁸
Tinnitus therapies	●	-	-	-	-
CROS/BiCROS	● ^{2,21}	-	-	-	-
Push button	-	● ²⁴	-	●	●
Volume control	● ²⁴	-	● ²⁴	● ²⁴	● ²⁴
Rechargeable	-	-	-	-	-
T-coil	○	-	●	●	●
Battery door on/off function	●	●	●	●	●
Audio input	● ²⁸	-	○ ³³	○ ³³	○ ³³
Modifications					
Color conversion kit	○	○	○	○	-
App/Accessory compatibility					
Signia app	●	●	●	●	●
Smart Connect	-	-	-	-	-
StreamLine TV/StreamLine Mic	●	-	-	-	-
miniPocket	●	●	●	●	●
easyPocket	-	-	-	-	-
ConnexxAir	-	-	-	-	-
Noahlink Wireless	●	-	-	-	-
Programming adapter/cable	size 13	size 312	size 13	Programming shoe/ CS 44 cable	Programming shoe/ CS 44 cable
Fitting set	ThinTube/Earhook	LifeTube	LifeTube	-	-

(2) Not available with PL 2/1
 (13) Ultra HD e2e
 (16) To be activated via rocker switch, Signia app or miniPocket
 (18) Basic remote tuning via 4-band-equalizer
 (19) Full live remote tuning via Connexx
 (24) Programmable via Connexx
 (28) Via direct streaming
 (33) With battery compartment modification

● Standard
 ○ Optional

Detailed feature overview BTEs

Type	Prompt S	Prompt P	Prompt SP	Run S
Performance level	Basic			
Picture				
Product details				
Matrix (2 ccm)	ThinTube/ LifeTube Earhook			
	45/124 dB -	- 70/134 dB	- 80/140 dB	45/124 dB -
Battery size	312	13	675	312
Power consumption	0.9 mA	1.4 mA	2.4 mA	0.9 mA
Features				
Frequency channels/ Compression controls	8/4	8/4	8/4	8
Hearing programs	4	4	4	6 preconfigured sound profiles
e2e wireless	-	-	-	-
Direct streaming	-	-	-	-
Push button	●	●	●	-
Volume control	-	● ⁽²⁴⁾	● ⁽²⁴⁾	●
T-coil	-	●	●	-
Battery door on/off function	●	●	●	●
Audio input	-	○ ⁽³³⁾	○ ⁽³³⁾	-
App/Accessory compatibility				
Signia app	●	●	●	●
Smart Connect	-	-	-	-
miniPocket	●	●	●	●
Programming adapter/cable	size 312	Programming shoe/ CS 44 cable	Programming shoe/ CS 44 cable	MobileFitting app
Fitting set	LifeTube	Comfort acoustic coupling	Comfort acoustic coupling	Comfort acoustic coupling

(24) Programmable via Connexx
(33) With battery compartment modification

● Standard
○ Optional

Run P	Run SP	Fun P	Fun SP	Fast P
Basic				
				
-	-	-	-	-
70/132 dB	79/140 dB	70/132 dB	79/140 dB	62/130 dB
13	675	13	675	13
1.3 mA	2.2 mA	1.2 mA	2.2 mA	1.1 mA
8	8	6	6	4
6 preconfigured sound profiles	6 preconfigured sound profiles	4 preconfigured sound profiles	4 preconfigured sound profiles	3 preconfigured sound profiles
-	-	-	-	-
-	-	-	-	-
-	-	-	-	-
●	●	●	●	●
-	-	-	-	-
●	●	●	●	●
-	-	-	-	-
●	●	-	-	-
-	-	-	-	-
●	●	-	-	-
MobileFitting app	MobileFitting app	Fitting via rocker switch	Fitting via rocker switch	Fitting via rocker switch
Comfort acoustic coupling	Comfort acoustic coupling	Comfort acoustic coupling	Comfort acoustic coupling	Simple acoustic coupling

Detailed feature overview Click ITEs and ITEs

Type	Active Pro	Active	Silk X Click CIC	Silk Nx Click CIC	Insio Nx IIC/CIC	Insio Nx ITC	Insio Nx ITE	Intuis 3 Click CIC	Intuis 3 Click ITC	Intuis 3 CIC	Intuis 3 ITC	Intuis 3 ITE
Performance levels			7X/5X/3X/2X/1X		7Nx/5Nx/3Nx/2Nx/1Nx			Essential		7px/5px/3px/2px/1px		
Picture												
Product details												
Matrix (2 ccm)	42/118 dB	42/118 dB	50/114 dB	50/113 dB	S 50/113 dB1 M 55/118 dB1 P 65/124 dB1	S 50/113 dB1 M 55/118 dB1 P 65/124 dB1	M 55/118 dB1 P 65/124 dB1	50/113 dB	55/116 dB	50/113 dB 55/118 dB 65/124 dB1	50/113 dB 55/118 dB 65/124 dB1	55/118 dB 65/124 dB1
Battery size	Lithium-ion power cell	Lithium-ion power cell	10	10	10	312	312	10	312	10	312	13
Power consumption			1.2 mA	1.3 mA	1.3 mA1	1.3 mA1	1.3 mA1	1.1 mA	1.1 mA	1.1 mA1	1.1 mA1	1.1 mA1
Features												
Frequency channels/ Compression controls	48/20	16/8	7X – 48/20 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	12/6	12/6	12/6	12/6	12/6
Hearing programs	6	6	6	6	6	6	6	4	4	4	4	4
e2e wireless	● ¹³	● ¹³	● ¹³	● ¹³	○ ¹³	○ ¹³	○ ¹³	-	-	-	-	-
Direct streaming	●	●	-	-	-	● ³⁴	● ³⁴	-	-	-	-	-
TwinPhone	● ¹⁵	-	● ¹⁵	● ¹⁵	● ³⁵	● ³⁵	● ³⁵	-	-	-	-	-
TeleCare	● ¹⁹	● ¹⁹	●	● ¹⁸	● ¹⁸	● ^{19, 34}	● ^{19, 34}	● ¹⁸				
Tinnitus therapies	●	-	●	●	●	●	●	-	-	-	-	-
CROS/BiCROS	-	-	●	● ³⁶	● ³⁶	● ²¹	● ²¹	-	-	-	-	-
Push button	○ ²⁵	○ ²⁵	○ ²⁵	-	○ ³⁷	○ ³⁴	○ ³⁴	-	●	○	○ ³⁴	○ ³⁴
Volume control	○ ²⁵	○ ²⁵	○ ²⁵	○ ²⁵	○ ²⁵	○ ²⁵	○ ²⁵	○ ²⁵	○ ²⁵	○ ²⁵	○ ²⁵	○ ²⁵
T-coil	-	-	-	-	-	-	-	-	-	-	○	○
Battery door on/off function	-	-	●	●	●	●	●	●	●	●	●	●
Audio input	● ²⁸	● ²⁸	-	-	-	○ ²⁸	○ ²⁸	-	-	-	-	-
Optivent	-	-	-	-	○	○	○	-	-	○	○	○
App/Accessory compatibility												
Signia app	●	●	●	●	●	●	●	●	●	●	●	●
Smart Connect	-	-	-	-	-	-	-	-	-	-	-	-
StreamLine TV/StreamLine Mic	●	●	-	-	-	● ³⁴	● ³⁴	-	-	-	-	-
miniPocket	●	●	●	●	●	●	●	●	●	●	●	●
easyPocket	-	-	-	-	-	-	-	-	-	-	-	-
ConnexxAir	-	-	●	● ³⁸	● ³⁸	-	-	-	-	-	-	-
Noahlink Wireless	●	●	-	-	-	●	●	-	-	-	-	-
Programming adapter/cable			Flex connector	Flex connector	Flex connector 37	Flex connector						

(1) Max. values
 (13) Ultra HD e2e with binaural audio signal processing
 (15) To be activated via Signia app or miniPocket
 (18) Basic remote tuning via 4-band-equalizer
 (19) Full live remote tuning via ConnexxAir
 (20) Not for 1px
 (21) With CROS Pure 312 Nx
 (25) Standard on Bluetooth-enabled ITC/ITE using Signia app/
 miniPocket, optional on other ITC/ITE

(28) Via direct streaming
 (29) With Smart Connect
 (34) Standard on Bluetooth-enabled ITC/ITE,
 optional on other ITC/ITE
 (35) To be activated via Signia app, miniPocket
 or optional push button
 (37) Not available for IIC
 (38) Only if e2e wireless is available

● Standard
 ○ Optional

Detailed feature overview Click ITEs and ITEs

Type	Prompt Click CIC	Prompt Click ITC	Prompt CIC	Prompt ITC	Prompt ITE	Run Click CIC	Run Click ITC
Performance levels	Basic					Basic	
Picture							
Product details							
Matrix (2 ccm)	50/113 dB	55/116 dB	50/113 dB 55/118 dB 65/124 dB 1	50/113 dB 55/118 dB 65/124 dB 1	55/118 dB 1 65/124 dB 1	50/113 dB	55/116 dB
Battery size	10	312	10	312/10	13/312	10	312
Power consumption	1.1 mA	1.1 mA	1.1 mA1	1.1 mA1	1.1 mA1	1.1 mA	1.1 mA
Features							
Frequency channels/ Compression controls	8/4	8/4	8/4	8/4	8/4	8	8
Hearing programs	4	4	4	4	4	6 preconfigured sound profiles	6 preconfigured sound profiles
Push button	-	●	○	○	○	-	-
Volume control	-	-	○ ²⁵	○ ²⁵	○ ²⁵	-	-
T-coil	-	-	-	○	○	-	-
Battery door on/off function	●	●	●	●	●	●	●
Optivent	-	-	○	○	○	-	-
App/Accessory compatibility							
Signia app	●	●	●	●	●	●	●
miniPocket	●	●	●	●	●	●	●
Programming adapter/cable	Flex connector	Flex connector	Flex connector	Flex connector	Flex connector	MobileFitting app	MobileFitting app

(1) Max. values
 (25) Standard on Bluetooth-enabled ITC/ITE using Signia app/
 miniPocket, optional on other ITC/ITE

● Standard
 ○ Optional