

Benefits Overview



Digital Assistant



Li-ion Rechargeability



Bluetooth Connectivity

Audiologist Benefits

- Combines advanced Signia Xperience audiology with the modern hearing design to combat stigma.
- Supports a wide range of patient needs, whether they're looking for all-day support or an occasional hearing boost.
- Modern earbud design that combats stigma and drives patient satisfaction with convenient, consumer-friendly features like Bluetooth® streaming, portable rechargeability and wireless charging.
- A great alternative to customs with Signia's Click Sleeves 2.0 – provides fast-fit technology in different sizes for same-day fitting and hearing aids that don't slip.

Key Features

- 4 hours for a full charge with up to 26 hours of wear time
- 30 minute quick charge for an average of 5 hours of wear time
- Automatic on/off functionality with charger
- Ultra HD e2e
- Remote Care Ready

Compatible with:



StreamLine TV™



StreamLine Mic™



miniPocket™

Wearer Benefits

- First-of-its-kind sleek, comfortable, and modern hearing aid design with the benefits of high-quality hearing without the stigma of a hearing aid.
- Pop them in when you need a hearing boost or keep them in all day long for non-stop hearing performance.
- Flexible charging options keeps you powered on the go – the portable charging case holds 3 days of extra charges and the built-in qi wireless charging technology provides added convenience.



Tinnitus



Recharge on-the-go



Signia App



Automatic Situation Detection



Brilliant Sound



Instant Fit



Enhanced Speech Understanding



Dynamic Soundscape Processing



White/Rose Gold



Black/Silver



Black

Signia App equipped with Signia Face Mask Mode

Improves audibility for those communicating with someone wearing a face mask - all with the touch of a button in the Signia App.



Type	Active Pro X
Matrix	42/118 dB
Battery Size	Lithium-ion power cell
Features	
Frequency Channels/Compression Controls	48/20
Hearing Programs	6
Bluetooth® Streaming	•
TwinPhone ¹	•
TeleCare	•
Tinnitus Therapies	•
CROS/BiCROS	--
Portable Rechargeability	•

Type	Active Pro X
T-coil	--
On/Off Function	Via Charger
App/Accessory Compatibility	
Signia App	•
Streamline TV/Streamline Mic	•
miniPocket	•
Noahlink Wireless	•
Programming adapter/cable line	--
Fitting set	Click Sleeves 2.0
Charging Options	USB-C and Qi wireless

(1) Activated Via Signia App or miniPocket.

The 'Qi' symbol is a trademark of the Wireless Power Consortium. The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners. SIV-10230 SH202115 Copyright © 2021 Signia GmbH. All rights reserved.

**Be
Brilliant™**