

Signia Motion Charge&Go X – Standard, Power and Super Power

# Benefits Overview



24/7 Digital Assistant



Bluetooth Connectivity



Li-ion Rechargeability

## Audiologist Benefits

- Complete lineup of rechargeable, Signia Xperience BTEs with the **world's first rechargeable Super Power** for added patient convenience.
- **Smaller, sleek designs** from previous models for added comfort.
- **Advanced benefits** with Power and Super Power options.
  - Dry&Clean Charger to **keep hearing aids in top working conditions** and reduce appointments.
  - **T-Coil** for patients leveraging looping systems for enhanced listening.

- **Compact chargers with long-lasting wear time** to get up to 60+ hours of wear time on a single charge. Charger models and options vary by model type.

## Key Features

- Up to 61 hours of use on a single charge with SP, 30 hours for P and 24 hours for Standard (based on a full 4-hour charge)
- 30 minute quick charge for an average of 5 hours of use
- Automatic on/off functionality with charger
- Extended frequency range offers wearers a competitive edge over other BTEs
- *Bluetooth®* connectivity
- Remote Care Ready via TeleCare 3.0

## Compatible with:



StreamLine TV™



StreamLine Mic™



miniPocket™

## Wearer Benefits

- Comfort and convenience without limits - **technology to support moderate, severe and profound hearing losses.**
- **Enhanced situational awareness** with adaptive technology that automatically adjusts to movements and changes around you.



Own Voice Processing



Dynamic Soundscape Processing



Acoustic Motion Sensors



Signia App




Tinnitus



T-coil

## Available Colors:



## Signia App is equipped with Signia Face

**Mask Mode** Improves audibility for those communicating with someone wearing a face mask - all with the touch of a button in the Signia App.

Type	Motion Charge&Go X	Motion Charge&Go Power X	Motion Charge&Go Super Power X
Available Performance Level	7	7	7
Matrix	60 dB/133 dB SPL – Standard Earhook 60 dB/125 dB SPL – Thin Tube 3.0	77 dB/135 dB SPL – Standard Earhook 66 dB/130 dB SPL – Thin Tube 3.0	82 dB/140 dB SPL
Battery Size	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell
<b>Features</b>			
Frequency channels/ compression controls	7X – 48/20	7X – 48/20	7X – 34/16
Hearing Programs	6	6	6
Bluetooth® Streaming	•	•	•
TwinPhone <sup>1</sup>	•	•	•
TeleCare 3.0	•	•	•
Tinnitus therapies	•	•	•
CROS/BiCROS <sup>2</sup>	•	•	•
Rocker Switch <sup>3</sup>	•	•	•
T-coil	--	•	•

(1) Activated via rocker switch, app or accessory

(2) Compatible with CROS Pure Charge&Go X

(3) Programmable via Connexx

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners. SIV-10231 SH202113 Copyright © 2021 Signia GmbH. All rights reserved.

Type	Motion Charge&Go X	Motion Charge&Go Power X	Motion Charge&Go Super Power X
<b>App/Accessory Compatibility</b>			
Signia App	•	•	•
StreamLine TV/ StreamLine Mic	•	•	•
miniPocket	•	•	•
Noahlink Wireless	•	•	•
Fitting set	EarWear 3.0	EarWear 3.0	Standard Tubing and Earmolds

All Motion Charge&Go X models work exclusively with the new EarWear 3.0. Contact your Signia Representative to learn more about the new EarWear 3.0.



Signia Motion Charge&Go X



Signia Dry&Clean Charger

One version compatible with the Power model, and a second option compatible with the Super Power model

**Be  
Brilliant™**