

Signia RIC product style and color guide

Pure® Charge&Go X



Made for
iPhone | iPad | iPod



Rechargeable

Pure® 312T X



312 Battery Type



Made for
iPhone | iPad | iPod



1	2	3	5	926	7	31	37	72	927
Beige	Granite	Grey	Black	Silver	Pearl White	Sandy Brown	Dark Champagne	Deep Brown	Rose Gold

Styletto™ X



Made for
iPhone | iPad | iPod



Rechargeable





				
Black/ Silver	Cosmic Blue/ Rose Gold	Snow White/ Rose Gold	Black	White

Actual colors may vary.

Click Sleeve 2.0


Products	Content	Item No.
Type	Size	
	xS	10951152
	S	10951146
	M	10951148
	L	10951150
	xS	10951153
	S	10951147
	M	10951149
	L	10951151

Click Domes


Products	Content	Item No.
Type	Size	
	4 mm	10426016
	8 mm	10426014
	10 mm	10426015
	8-12 mm	10670790
	6 mm	10426017
	8 mm	10426012
	10 mm	10426013
	8/10 mm	10426018
	10/12 mm	10426019

Fitting Ranges


S-Receiver Fitting Range 45/108
– 45 dB Receiver




M-Receiver Fitting Range
60/119 – 60 dB Receiver











P-Receiver Fitting Range
70/124 – 70 dB Receiver



Integrated HP-Receiver Fitting Range
75/130 – 75 dB Receiver
(Not available for Styletto X)



Common RIC Earmolds

									
	ROES Codes	RIC	Embedded RIC¹	RIC Lock	Embedded RIC Lock¹	RIC Foil²	Foil RIC with Lock²	Canal RIC	Canal Lock RIC
		EM11	EM28	EM12	EM29	EM15	EM16	EM17	EM18
Material	Acrylic (Hard)	•	•	•	•	•	•	•	•
	Silicone (Soft)	•	•	•	•	—	—	•	•
Vent Type	No Vent	•	•	•	•	•	•	•	•
	Semi-IROS	•	•	•	•	•	•	•	•
	Standard	•	•	•	•	•	•	•	•
Vent Size	1.0mm	•	•	•	•	•	•	•	•
	1.2mm	•	•	•	•	•	•	•	•
	1.6mm	•	•	•	•	•	•	•	•
	2.0mm	•	•	•	•	•	•	•	•
	3.0mm	•	•	•	•	•	•	•	•
	Largest Vent Possible	•	•	•	•	•	•	•	•

1) Embedded fixed receiver 2) Acrylic material only, applies to the Gloss

Common Earmold Color Options

Gloss






Matte



Signia App



-  Signia Assistant
-  Face Mask Mode
-  TeleCare



Accessories



StreamLine TV™



StreamLine Mic™



miniPocket™

The **Bluetooth®** word mark and logos are owned by the **Bluetooth SIG, Inc.** and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners.

Android is a trademark of Google Inc.

iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.

"Made for iPhone" means that an electronic accessory has been designed to connect specifically to iPhone and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards.

Please note that the use of this accessory with iPhone may affect wireless performance.