

# Signia Pocket Guide

**Fall 2021** – The most comprehensive rechargeable portfolio with  
the only contactless custom rechargeable hearing aid

[veterans.signiausa.com](https://veterans.signiausa.com)

## Signia Augmented Xperience Platform and Signia Xperience Platform



Insio Charge&Go AX  
(Custom)



Pure® Charge&Go T AX  
(RIC)



Stylecto X  
(RIC)



Pure® 312T X  
(RIC)



Motion Charge&Go X (BTE)  
Standard, Power and Super Power

## Expanded Portfolio Offerings



NEW

**Insio  
Charge&Go  
AX**



NEW

**Pure  
Charge&Go  
T AX**



NEW

**CROS Pure  
Charge&Go  
AX**



**Styletto X**



**Pure 312T  
X**



**Motion  
Charge&Go  
X**



**Motion  
Charge&Go  
PX**



**Motion  
Charge&Go  
SPX**

Signia is dedicated to improving the quality of life for all veterans, active military, and American Indians through better hearing. Through innovative solutions, unparalleled support, and service excellence, Signia is proud to work with Government Services audiologists to enhance their patients' human performance and help them experience the world around them on their terms."

The Signia AX platform is a ground-breaking audiological advancement and the world's first to use two separate processors that automatically work together to deliver outstanding speech clarity in quiet, moderate, or noisy environments. In studies, patients describe 25 percent greater speech understanding in noise when they use AX-based hearing aids; 95 percent of study participants reported exceptional speech understanding in their home environment.

In addition to these new products, Signia reaffirms its commitment to Government Services audiologists and their patients through its Signia Delivers promise. Sharing their goal of achieving the best patient outcomes possible, Signia Delivers across three pillars:

- **Innovative Solutions:** In addition to its robust product offerings, including the world's largest portfolio of rechargeable solutions, Signia has introduced pioneering solutions like Signia Remote Care for virtual hearing care and adjustments and Face Mask Mode.
- **Unparalleled Support:** the company's recently expanded team of audiology, customer care, and tech support professionals brings an unrivaled 180 combined years of Government Services expertise to understand audiologists' needs and assist with product selection, challenging fittings, personalized training, and more.
- **Service Excellence:** Signia strives to handle each order as efficiently and professionally as possible, maintaining a laser focus on enhancing order fulfillment and reducing delivery time.

# Table of Contents

## Our People.....5

### Audiology

Signia Augmented Xperience .....	11
Signia Xperience.....	14
Signia App.....	15
Signia TeleCare.....	16
Connexx.....	17
Signia Fitting .....	18

### Technology

Insio Charge&Go AX.....	21
Pure Charge&Go T AX.....	23
CROS RICs.....	24
Signia AX RIC Fitting Ranges.....	25
Receiver 3.0.....	26
EarWear 3.0.....	27
Styletto X.....	28
Signia X RIC Fitting Ranges.....	29
Motion Charge&Go X.....	30
Motion Charge&Go X fitting ranges .....	31
Thin Tube 3.0 .....	32
EarWear 3.0.....	33
Pure 312 X.....	34
Insio Nx Customs .....	35
Insio Nx Custom products and fitting ranges.....	36
Motion 13 and 13 P Nx.....	37

### Streamers

StreamLine™ Mic.....	39
StreamLine™ TV.....	42

### Remote Controls

miniPocket™ .....	45
Pairing hearing aids with miniPocket.....	46

### Chargers

Charger Compatibility .....	47
Rechargeability.....	49
Portable AX Charger.....	50
Insio Charge&Go AX Charger.....	51
Closed-lid Inductive Charger.....	52
Styletto X Charger.....	53
Dry&Clean Charger.....	54
How to order a Charger in ROES.....	56

### Earmolds

RIC Earmolds.....	57
BTE Earmolds .....	58

## ROES Codes .....59





## Government Services Team

**Phone:** (800) 955-4327 • **Fax:** (732) 562-6653

**Consumer Hotline:** (800) 350-6093

**WEST**

Alaska  
Hawaii

WA, OR, ID, MT, ND, WY, SD, NE, UT, CO, NM, AZ, CA, NV, TX-West, TX-East

**CENTRAL**

MN, WI, MI, IA, IL, IN, OH, KY, PA, NY, VT, NH, ME, CT, RI, MA, NJ, DE, MD, WV, VA, NC, SC, GA, FL

**SOUTH**

LA, MS, AL, TN, MO, AR, OK, KS

**MOUNTAIN**

**SOUTHEAST**

**Puerto Rico**

**Team Members:**

**WEST:** Jessica Hillham, Ingrid McBride, Daniela Valencia

**CENTRAL:** Peggy Zannmiller, Jennifer Scheytt, Kenya Barley

**SOUTH:** Keri Morgan, JoAnn Kerwin, Kenya Barley

**MOUNTAIN:** Daniel Greenblatt, Crystal Dvorack, Daniela Valencia

**SOUTHEAST:** Spencer Reeder, Courtney Rogers, Sarah Andresen

**NORTHEAST:** Jim Coyne, Rachel Pokorney, Sarah Andresen

## Your Signia Team

**Phone:** (800) 955-4327 • **Fax:** (732) 562-6653

**Consumer Hotline:** (800) 350-6093

### Northeast Team



**Government Account Manager**

Jim Coyne, Au.D.

(585) 236-0636

james.coyne@signiausa.com



**Government Audiologist**

Rachel Pokorney, Au.D.

rachel.pokorney@signiausa.com



**Customer Care**

Sarah Andresen

sarah.andresen@signiausa.com

### Southeast Team



**Government Account Manager**

Spencer Reeder, Au.D.

spencer.reeder@signiausa.com



**Government Audiologist**

Courtney Rogers, Au.D.

courtney.rogers@signiausa.com



**Customer Care**

Sarah Andresen

sarah.andresen@signiausa.com

## Your Signia Team

**Phone:** (800) 955-4327 • **Fax:** (732) 562-6653

**Consumer Hotline:** (800) 350-6093

### South Team



**Government Account Manager**

Keri Morgan, Au.D.

(417) 291-1411

[keri.morgan@signiausa.com](mailto:keri.morgan@signiausa.com)



**Government Audiologist**

JoAnn Kerwin, Au.D.

[joann.kerwin@signiausa.com](mailto:joann.kerwin@signiausa.com)



**Customer Care**

Kenya Barley

[kenya.barley@signiausa.com](mailto:kenya.barley@signiausa.com)

### Central Team



**Government Account Manager**

Peggy Zanmiller, B.S. Speech and Hearing Science

(952) 221-8379

[margaret.zanmiller@signiausa.com](mailto:margaret.zanmiller@signiausa.com)



**Government Audiologist**

Jennifer Scheytt, Au.D.

[jennifer.scheytt@signiausa.com](mailto:jennifer.scheytt@signiausa.com)



**Customer Care**

Kenya Barley

[kenya.barley@signiausa.com](mailto:kenya.barley@signiausa.com)

## Your Signia Team

**Phone:** (800) 955-4327 • **Fax:** (732) 562-6653

**Consumer Hotline:** (800) 350-6093

### Mountain Team



**Government Account Manager**

Daniel Greenblatt, Au.D.

(480) 249-8974


[daniel.greenblatt@signiausa.com](mailto:daniel.greenblatt@signiausa.com)



**Government Audiologist**

Crystal Dvorak, Au.D.

[crystal.dvorak@signiausa.com](mailto:crystal.dvorak@signiausa.com)



**Customer Care**

Daniela Valencia

[daniela.valencia@signiausa.com](mailto:daniela.valencia@signiausa.com)

### West Team



**Government Account Manager**

Jessica Hillham, Au.D.

[jessica.hillham@signiausa.com](mailto:jessica.hillham@signiausa.com)



**Government Audiologist**

Ingrid McBride, Au.D.

[ingrid.mcbride@signiausa.com](mailto:ingrid.mcbride@signiausa.com)



**Customer Care**

Daniela Valencia

[daniela.valencia@signiausa.com](mailto:daniela.valencia@signiausa.com)

## Your Signia Team

**Phone:** (800) 955-4327 • **Fax:** (732) 562-6653

**Consumer Hotline:** (800) 350-6093



**VP Government Services**

John Murray  
[john.murray@signiausa.com](mailto:john.murray@signiausa.com)



**Lead Engineer**

Andrei Bucur  
[andrei.bucur@signiausa.com](mailto:andrei.bucur@signiausa.com)



**Director of Sales**

Brian Fasler  
[brian.fasler@signiausa.com](mailto:brian.fasler@signiausa.com)



**Director of Product Management**

Julie Dunphy  
[julie.dunphy@wsa.com](mailto:julie.dunphy@wsa.com)



**Manager of Government Services**

Kerrin McDevitt  
[kerrin.mcdevitt@signiausa.com](mailto:kerrin.mcdevitt@signiausa.com)



**Contract Manager**

Yazmin Mendoza  
[yazmin.mendoza@signiausa.com](mailto:yazmin.mendoza@signiausa.com)







## Augmented Xperience Platform at a Glance



2x

Audiological core  
processors

110%

Increased binaural  
audio transfer rate

60%

Reduced binaural  
latency

27%

More electric  
circuits

22%

Longer battery  
runtime

## Augmented Xperience Platform at a Glance

### New Audiological Features



#### e2e wireless 4.0: A great leap for true binaural hearing

e2e 4.0 significantly boosts the speed of data transmission and a 60% reduction in latency provides inaudible shifts between environments. The synchronization of more binaural audio data than ever before ensures clearer speech processing with high-precision directional awareness.



#### Increased input dynamic range

The microphones and the state-of-the-art analogue input stage can handle signals up to 117dB with great sound quality. That's very close to the full capacity of human ears and ensures a clean, crisp signal for scenarios like concerts.



#### TruEar 360

The powerful processing engine allows for an exceptional TruEar 360 mic mode that reproduces and emphasizes the capabilities of the natural pinna-effect better than ever. TruEar 360 highlights the difference between front and back sources of sound to make localization of sound sources easier.



#### Enhanced feedback cancellation

e2e 4.0 also synchronizes feedback cancellation between the ears in the same phase, removing more artifacts for a smoother overall sound experience.

## Augmented Xperience Platform at a Glance

### Enhanced Audiological Features



#### Direct streaming for iOS and Android

Augmented Xperience offers direct streaming for both iOS and ASHA-compatible Android smartphones. It handles the signal of a streaming input independently from other signals, ensuring that the sound is even clearer and easier to optimize.



#### Automatic Situation Detection

Provides detailed information about the wearer's movement and environment. Analyzing and adapting to the acoustic and motion changes, these sensors deliver a clear understanding of speech from any direction and maximizing audibility as you move throughout the day.



#### Own Voice Processing (OVP™)

Can separate the wearer's own voice from other speech and noise for a natural sound, increasing first-fit acceptance to 80% and keeping patients engaged in conversations.<sup>1</sup>

\*Excluding Insio Charge&Go

## Signia Xperience Platform

Signia Xperience is the industry's leading audiological platform, and it powers Signia hearing aids. It focuses on the lifestyle and actions of the wearer for optimized support in every situation. Intuitive and flexible programming software makes Signia Xperience adaptable and user-friendly, just another way we're putting the wearer first.

AMS

### Acoustic Sensors

Analyzing and adapting to the acoustic changes around you, these sensors deliver personalized amplification and a clear understanding of speech from any direction, mirroring optimized natural hearing experiences.

### Motion Sensors

Optimizing hearing in every situation by recognizing changes in your own movement to adjust acoustic settings and maximizing audibility as you move throughout the day.

OVP

### Own Voice Processing (OVP™)

Separates the wearers own voice from other speech and noise for a natural sound, increasing first-fit acceptance to 80% and keeping patients engaged in conversations.

DSP

### Dynamic Soundscape Processing

Utilizing Signia's iconic binaural audio processing with Binaural Directionality, this world's first technology has been clinically proven to provide better than normal hearing in noisy situations.

1 Haubold J. (2019) Real-world effectiveness of Signia Xperience Dynamic Soundscape Processing. Unpublished Manuscript.

2 Froehlich M, Branda E, Freels K. (2019) New dimensions in automatic steering for hearing aids: Clinical and real-world findings. Hearing Review. 26(11):32-36.

3 Høydal, E. H. (2017). A new own voice processing system for optimizing communication. Hearing Review, 24(11), 20-22. Further details: [www.signia-pro.com/ovp-study](http://www.signia-pro.com/ovp-study).

4 Froehlich M, Freels K, Powers T. (2015). Speech recognition benefit obtained from binaural beamforming hearing aids: comparison to omnidirectional and individuals with normal hearing AudiologyOnline, Article 14338. Retrieved from <http://www.audiologyonline.com>

## Signia App for Patients

Beyond streaming and remote control capabilities, the Signia app is customized to display specific features based on the patient's hearing aid model and gives them freedom and support with access to:




**Telecare** for remote care support with you! Patients fitted with Signia Telecare have a 19% higher acceptance rate than patients fitted without it.



**Signia Assistant** for on-demand AI support. 93% stated that using the Signia Assistant improved their satisfaction for listening in difficult situations.



**Signia Face Mask Mode** for an immediate boost in speech with the simple tap of a button. Facemasks have been shown to reduce audibility by up to 30%, and make lip reading impossible.



## Signia TeleCare – 3.0

Proven remote care support to remain virtually connected to your patients

**Remain connected with patients — Communicate with patients through Signia's safe and secure connection on any laptop, smartphone, or computer.**

**Higher patient satisfaction — More flexibility to support patient needs.**

- In-person or remote activation of TeleCare for new and existing patients.
- Remote adjustments via Connexx.
- Compatible with multiple Signia platforms (Nx, X).
- Adjustments are sent to the patient in real-time (synchronously).



Patients fitted with Signia Telecare have a 19% higher acceptance rate than patients fitted without it.\*



\* Article, Matthias Froehlich, PhD, Eric Branda, AuD, PhD, Dirk Apel, BS, 2020: User Engagement with Signia TeleCare: A Way to Facilitate Hearing Aid Acceptance

## Connexx 9.6.0 UPDATE

- Improved performance speed of software.
- Wireless hearing aid firmware updates are available via ConnexxAir wireless, only for non-Bluetooth® hearing aids. Firmware updates are only available when received from the notification center.
- And now Insitu-Audiometry available through remote TeleCare.



version **9.6.0** with SIFIT™ 9.6.0

# Signia Fitting Guide

**NOTE: All ordered hearing aids have already been First Fit.**

1. Select your patient from NOAH patient database
2. Launch Signia Connexx fitting software
3. Select "Detect Hearing aids"



Audiogram




Detect Hearing Instruments



Hearing Loss Simulator

4. Once hearing aids are detected, read out connected instruments




5. Run OVP (Own Voice Processing)

For optimal results



Complete the OVP training: Tell your patient that they will hear a beep, and that they should count loudly and continuously from 2100 upwards. Example 2100, 2200, 2300 and so on.



OVP is automatically enabled. Select "OK".  Training was successful

Continued on next page...

## Signia Quick Fit Guide

(continued)

6. Make any gain adjustments with Basic or Fine tuning.

**Basic Tuning**

Basic Tuning

**Fine Tuning**

Fine Tuning

7. Select **Configuration**

Configuration

- **Hearing aid** – Configure User Controls, Volume Control, Therapy Signal, SoundBalance.
- **System Sounds** – Configure Alert Tones.

8. Save Session and Select "X" at top right of screen or select drop down item under File.

- Disconnect hearing aids.

## CROS Pure Charge&Go X and CROS Pure 312T X Quick Fit Guide

Under Fine Tuning, select Sound Settings to adjust the sensitivity of the CROS microphone in relation to the hearing aid mic.

If the “Enable CROS microphone” is unchecked, the CROS microphone will be disabled.

### Input Configuration

Mode: CROS

☒ Enable CROS microphone

CROS level



Adjust the CROS microphone level





# Insio Charge&Go AX Customs



Note: if full shell model is needed, please order a half shell and select "full shell dexterity" under dexterity options



**Insio Charge&Go AX Charger**  
(Standard) – features contactless charging

## Key features & benefits



Rechargeable



iPhone &  
Android  
Connectivity



CROS



Tinnitus



Augmented  
Speech  
Understanding



Immersive  
Soundscape



Automatic  
Situation  
Detection



Signia App



AI Digital  
Assistant



Black



Dark Brown



Brown



Mocha






Tan



Beige

Please note that red, blue and clear shells are also available.

## Insio Charge&Go AX Customs



## Pure® Charge&Go T AX

The ultimate rechargeable RIC.



### Portable Charger

(Standard) - 36 hours  
on a single charge



Black



Graphite



Dark  
Champagne



Silver



Pearl  
White



Fine  
Gold



Deep  
Brown



Sandy  
Brown



Rose  
Gold



Beige

### Key features & benefits



Integrated  
Telecoil



Android  
& iPhone  
Connectivity



Recharge  
on-the-go



CROS



Augmented  
Speech  
Understanding



Immersive  
Soundscape



Automatic  
Situation  
Detection



AI Digital  
Assistant

- Receivers: 3.0 S, M, P, HP

- EarWear: EarWear 3.0

- Additional Chargers Available: Pure Charger,  
Dry&Clean Charger

## CROS — RICs



**CROS Pure  
Charge&Go T AX**



**CROS Pure  
Charge&Go X**



**CROS Pure  
312X**

- Li-ion rechargeability and auto on/auto off function when removed from charger for easy handling.
- World's first CROS/BiCROS technology to benefit from acoustic-motion sensors for the clearest speech understanding that adapts to the wearer's movements (with BiCROS).
- Own Voice Processing (OVP™) (with BiCROS).
- Directional pattern skewed toward the unaidable ear for improved speech understanding.
- CROS mixing-level adjustment available in the Signia App to balance spatial benefit vs. noise level.
- Compatible with Signia App and miniPocket™.

### DID YOU KNOW?

Signia leads the industry in CROS/BiCROS technology.

**A survey\* among hearing care professionals found that:**

- Signia CROS technology received an average of **4.5 out of 5 stars**, the majority stating it is the **BEST TECHNOLOGY** on the market.
- Most reported that they and their patients were **EXTREMELY SATISFIED** with Signia CROS.

**When asked for specific feedback, here's what your peers had to say:**


- "It is **the only CROS product that I have had success with**. These patients still have some challenges, but overall they **accept the fitting at a much higher rate than any other CROS product** I have tried."
- "Patient is **hearing in situations they never heard** with their old [competitive CROS product]."
- Although we have used [competitor] in the past, my **previous [competitor] users who now have the Signia CROS are happier**.

## Signia AX RIC products and fitting ranges


### Pure Charge&Go T AX




- Lithium-ion power cell



Eartip 3.0 Open  
 + Sleeve 3.0 Power  
 + Earmold 3.0




Eartip 3.0 Open  
 + Sleeve 3.0 Power  
 + Earmold 3.0



Sleeve 3.0 Power  
 + Earmold 3.0

**75/131 – HP-Receiver**  
(Not available with Styletto X)



Custom Shell (no vent)

## Receiver 3.0

receiver types

Our new EarWear 3.0 offers a higher physical fit rate, improved cosmetics and better retention



miniReceiver 3.0 Cable Bending  
10 / 10 fit correctly – out of the box

Only available for Pure Charge&Go T AX  
and CROS Pure Charge&Go T AX

Smaller plug for smaller RICs

Push-/pull system –  
no tool required

Smart Receiver Detection  
prevents misuse


6-Pin Connector for  
new features

New bending ensures  
higher retention in ear

New bending sits closer  
= cosmetically more  
appealing

Size 0-4

Reinforced cable outlet for  
highest durability





## **NEW** EarWear 3.0

This versatile new range is suitable for the vast majority of ear anatomies and hearing needs

Soft eartips allow for deep insertion with high comfort - also into small ear canals

- The new Power Sleeve 3.0 expands the acoustically highly reliable sleeve to also offer a fully closed option.
- Secure, proven design makes Eartips 3.0 and Sleeves 3.0 much easier for wearers to exchange.
- All Earwear 3.0 items can be used together.



## Styletto X

Exchangeable receivers, rechargeability, and *Bluetooth®* streaming in a one-of-a-kind design.



reddot design award  
winner 2019



**Portable Charger**  
(Standard) – Holds  
3 days of charging

Available colors



Black/  
Silver

Cosmic Blue/  
Rose Gold

Snow White/  
Rose Gold

Black

Snow  
White

### Qi Wireless Technology<sup>1</sup>

Charger equipped with Qi Wireless Technology — enables wireless charging on any standard charging pad, including those built into cars and cafe counters.

### Key features & benefits



Recharge  
& Go



iPhone &  
Android<sup>2</sup>  
Connectivity



Tinnitus



Dynamic  
Soundscape  
Processing



Acoustic  
Motion Sensors



Own Voice  
Processing



Signia App



Face Mask  
Mode



AI Digital  
Assistant

- Receivers: miniReceivers 2.0 S, M, P

<sup>1</sup> Wireless charging pad not included with purchase

<sup>2</sup> Requires StreamLine Mic

## Signia X RIC products and fitting ranges

### Pure 312T X




- Battery Type 312

### Styletto X



- Lithium-ion power cell


#### 46/110 – S-Receiver



Open Click Domes  
Double Click Domes  
Click Mold (no vent)

Channels  
**48**  
Programs  
**6**


#### 60/119 – M-Receiver



Open Click Domes  
Double Click Domes  
Click Mold (no vent)

Channels  
**48**  
Programs  
**6**


#### 65/122 – P-Receiver



Double Click Domes  
Click Mold (no vent)

Channels  
**48**  
Programs  
**6**

#### 75/131 – HP-Receiver (Not available with Styletto X)



HP-Receiver

Channels  
**48**  
Programs  
**6**

## Motion Charge&Go X

Complete BTE lineup with the world's first rechargeable Super Power

New EarWear 3.0 thin tubes, domes and sleeves are required for thin tube fittings



Standard

Power

Super Power



Signia Motion Charge&Go X



Signia Dry&Clean Charger for P & SP

Up to 61 hours of use on a single charge with SP, 30 hours for P and 24 hours for Standard (based on a full 4-hour charge)

### Key features & benefits



Lithium-Ion



iPhone & Android<sup>1</sup> Connectivity



Integrated Telecoil<sup>1</sup>



CROS



Tinnitus



Dynamic Soundscape Processing



Acoustic Motion Sensors



Own Voice Processing



Signia App



Face Mask Mode




AI Digital Assistant

<sup>1</sup> P and SP models


Beige	Granite	Grey	Black	Silver Hero	Pearl White	Sandy Brown	Dark Champagne	Deep Brown	Rose Gold
1	2	3	5	926	7	31	37	72	927

## Motion Charge&Go X fitting ranges


Motion C&G X Fitting Range



Motion C&G P X Fitting Range



Motion C&G SP X Fitting Range





## Thin Tube 3.0

Compatible with Motion Charge&Go X and P BTE

Both ThinTube 3.0 options offer proven better retention, comfort, usability and physical fit rate

The new 1.4mm-width PowerTube also delivers a cosmetically appealing solution for more severe degrees of hearing loss.



### Retention

Improved bending of tube + outlet locks the acoustical coupling reliably into the ear canal, prevents lateral migration and makes concha-locks redundant.



Binaural Thin tube 3.0 S0	10991609
Binaural Thin tube 3.0 S1	10991610
Binaural Thin tube 3.0 S2	10991611
Binaural Thin tube 3.0 S3	10991612
Binaural Thin tube 3.0 S4	10991613
Binaural Thin tube 3.0 S0 Power	10992492
Binaural Thin tube 3.0 S1 Power	10992493
Binaural Thin tube 3.0 S2 Power	10992494
Binaural Thin tube 3.0 S3 Power	10992495
Binaural Thin tube 3.0 S4 Power	10992496

## **NEW** EarWear 3.0

This versatile new range is suitable for the vast majority of ear anatomies and hearing needs

Soft eartips allow for deep insertion with high comfort - also into small ear canals

- The new Power Sleeve 3.0 expands the acoustically highly reliable sleeve to also offer a fully closed option.
- Secure, proven design makes Eartips 3.0 and Sleeves 3.0 much easier for wearers to exchange.
- All Earwear 3.0 items can be used together.



## Pure® 312T X

Streaming RIC with T-Coil



### Key features & benefits



312 Battery



iPhone & Android<sup>1</sup> Connectivity



Integrated Telecoil



CROS



Tinnitus



Dynamic Soundscape Processing



Acoustic Motion Sensors



Own Voice Processing



Signia App



Face Mask Mode




AI Digital Assistant

- Receivers: miniReceivers 3.0 S, M, P, HP

Beige	Granite	Grey	Black	Silver Hero	Pearl White	Sandy Brown	Dark Champagne	Deep Brown	Rose Gold
1	2	3	5	926	7	31	37	72	927

<sup>1</sup> Requires StreamLine Mic

## Insio Nx Custom



### Custom Faceplate Colors NEW Non-Bluetooth® Faceplate available

Description	ITE	ITC/HS	CIC	IIC
Beige	•	•		
Tan	•	•		
Mocha	•	•	•	
Brown	•	•	•	
Dark Brown	•	•		
Black				•

### Custom Shell Colors

Description	ITE	ITC/HS	CIC	IIC
Beige	•	•		
Tan	•	•	•	
Brown	•	•	•	
Dark Brown	•	•		
Clear	•	•	•	
Red	•	•	•	•
Blue	•	•	•	•

### Programmers



Noahlink Wireless  
for ITC/HS and ITE





ConnexxAir  
for IIC and CIC



Stream audio directly into  
your *Bluetooth®* enabled Insio  
ITC/HS and ITE instruments  
with **Streamline™ TV** and  
**Streamline™ Mic**

## Insio Nx Custom



## Motion 13 and 13 P Nx



Motion 13P Nx



Motion 13 Nx



**Fitting Ranges**



Beige	Granite	Grey	Black	Silver	Pearl White	Golden Blonde	Sandy Brown	Dark Champagne	Dark Granite	Deep Brown
1	2	3	5	6	7	30	31	37	66	72








## StreamLine Mic – Compatible with all Signia Bluetooth® hearing aids

The new StreamLine™ Mic can be connected to virtually any Bluetooth®-enabled smartphone (such as iPhone and Android) and establishes various direct audio streaming options:

- Phone call streaming (both ears)
- Integrated Microphone
- Hands-free phone calls (both ears)
- Accept/decline calls
- Remote microphone functionality (both ears)
- Hear in challenging situations with remote microphone mode
- Audio, video and music streaming (both ears, stereo)
- Remote Control functions
- Volume Increase
- Volume Decrease
- Program Change

**Watch handling videos online**

[veterans.signiausa.com/howto](https://veterans.signiausa.com/howto)



### HOW TO ORDER IN ROES

**ROES Main Menu – Patient Orders.**

**Select an Action:** Wireless Devices

Remotes

Wireless Transmitter

**Select Device Type:** SIVANTOS

**Select Make:** StreamLine Mic

## StreamLine Mic – Overview

### Multi-function button (MFB)

- Switch on (press 3 sec) and off (press 6 sec)
- Toggle between two modes:
  - Bluetooth® Phone Mode
  - Remote Microphone Mode


### Wearer controls/status indicators

- VC / Program up
- VC / Program down
- Power indicator
- Status indicator






### StreamLine Mic offers intuitive onboard control and status elements

### Accessories included in the shipment

- 2 lanyards (regular and long)
- USB wall plug + USB connector
- Pouch
- User Guide



### StreamLine Mic controls

- 1 Microphones
- 2 Power indicator 
- 3 Status indicator 
- 4 Lapel clip
- 5 Multi-function button (volume, program, mute) 
- 6 Volume buttons  and 
- 7 USB port for charging

# StreamLine Mic – Pairing Quick Start Guide

## Pairing with Hearing Aids

- Open/close the battery doors of the hearing aids and/or place in/out of the charger to place them in pairing mode.
- Turn on StreamLine Mic by holding down the MFB until the power indicator is green. This takes approximately 3 seconds. The status indicator flashes yellow.
- Place both hearing aids close to StreamLine Mic (maximum distance 20 cm or 8 inches). The pairing process starts automatically. Pairing is finished when the status indicator is solid yellow.
- A solid green power indicator and solid yellow status indicator verify connection to hearing aids is successful with no mobile device paired.
- Your hearing aids need to be paired once with StreamLine Mic.

## Pairing with Smartphone

- Confirm that Bluetooth® of the mobile device is on.
- Turn on StreamLine Mic by holding down the MFB until the power indicator is green. This takes approximately 3 seconds.
- Hold down the the MFB and “+” button on the StreamLine Mic at the same time until status indicator flashes blue. This takes approximately 6 seconds.
- The wearer has 3 minutes to complete the pairing.
- Select “StreamLine Mic” from the list of detected devices on the Bluetooth® menu.
- If requested, enter the PIN “0000”.
- A solid green power indicator and solid blue status indicator verify connection to smartphone is successful.



## StreamLine Mic – Battery Runtime

- Charging status is clearly displayed via power indicator on the StreamLine Mic
- Full charge takes approx. 2 hours
- Up to 6 hours battery runtime after full charge with full functionality (phone calls, audio streaming, remote mic)
- Can be charged via USB or usual plug-in (all necessary cables are included in the delivery)
- Can be used while charging

## StreamLine TV – Compatible with all Signia Bluetooth® hearing aids

The perfect companion  
for TV fans.

### StreamLine TV.

Thanks to Streamline™ TV audio from television or any device with audio output can be fed directly into Signia X and Nx hearing aids. Wearers can even adjust the volume individually via the Signia app.

All wearers of AX, X and Nx hearing aids have direct binaural streaming for TV via the StreamLine TV.

This includes use without the Signia app for non-smartphone users.



#### HOW TO ORDER IN ROES

##### ROES Main Menu – Patient Orders.

**Select an Action:** Wireless Devices

Remotes

Wireless Transmitter

**Select Device Type:** SIVANTOS

**Select Make:** StreamLine TV

## Pairing for StreamLine TV

**If ordered together, Signia will pair all Bluetooth® compatible hearing aids to the Streamline TV.**

The hearing aids need to be paired once with the transmitter. After that, the pairing is saved in the hearing aids and remembered each time they are turned on.


- Ensure that the transmitter is connected to the power supply.
- Turn the hearing aids off and on again: open and close the battery compartments. If your hearing aids do not have a battery compartment, turn them off and on as described in the hearing aids' user guide.
- Place both hearing aids close to, or on top of, the transmitter. The pairing process starts automatically.
- Observe the indicators:



### Description



If pairing is completed, both indicators are blue for about 30 seconds.



### Multiple users:

- The transmitter can be paired to multiple hearing aids. All hearing aids will be able to hear the same audio signal.
- To clear the pairings of one current user, first reset the transmitter to factory settings. This will clear the pairings of all users. After that, pair the hearing aids of the other users again. Refer to the section "Resetting to factory settings".

## Connecting StreamLine TV

### TOSLINK connection

- Before connecting the cable, remove the protection caps from the TOSLINK cable plugs.
- Plug the cable into the TOSLINK **input** of the transmitter and into the TOSLINK (optical audio) **output** of electronic device. To prevent the TOSLINK cable from breaking, do not forcefully bend it.




### Further information

- Refer to the electronic device's manual: set the output format of the electronic device to either PCM (pulse-code modulation) or Dolby® Digital.\*

### RCA connection

- Plug the cable into the RCA connections of the transmitter.
- Plug the other side of the cable into the RCA output of the electronic device or via adapter into the SCART or via mini jack output of the electronic device.



For more information and details on how to connect to your home theater amplifier or a sound bar watch the handling videos online at [veterans.signiausa.com/howto](http://veterans.signiausa.com/howto)

\*Dolby and the double-D symbol are registered trademarks of Dolby Laboratories.





## miniPocket — Remote control compatible with all Signia hearing aids


Small enough to hang on a keyring, the miniPocket™ is ideal for wearers who want quick, discreet control of their hearing aids without using a smartphone.

- Battery type CR2450
- Pairing remains in place when battery dies and is replaced

---

### Components

- ❶ Status LED
  - ❷ Volume buttons
  - ❸ Keyring
  - ❹ Key lock switch
  - ❺ Program change button
- 



## Pairing hearing aids with miniPocket

1. Move key lock switch to 'lock' position



2. Turn on hearing aids and place in ears
3. Hold down **P** button on the remote control while switching the key lock to the unlock position
4. Release **P** button (LED on)
5. Press and hold **P** button until the pairing confirmation tone is heard
6. Move Key Lock switch to Lock position to exit pairing mode



**Watch video online**

[veterans.signiausa.com/howto](https://veterans.signiausa.com/howto)

### HOW TO ORDER IN ROES

When ordering a hearing aid in ROES, look for the "Select a remote" dropdown. Select "miniPocket"






## Charger Compatibility

Chargers							
	Portable Pure Charge&Go AX Charger	Pure Charge&Go AX Charger	Dry&Clean Pure Charge&Go AX Charger	Inductive Charger I	Inductive Charger II	Styletto X Charger	Styletto Connect Charger
type of charger	Galvanic	Galvanic	Galvanic	Inductive	Inductive	Galvanic	Galvanic
part number	10993118	10991650	10996345	10948199	10966942	10972015	10965237
      							
cord connection	USB-C	Attached	USB-C	micro USB	micro USB	USB-C	micro USB
Hearing Aids							
RIC	Pure C&G AX	* *	•	•			
	CROS Pure C&G AX	* *	•	•			
	Pure C&G X			•	* *		
	CROS Pure C&G X			•	* *		
	Styletto X					* *	
	Pure C&G Nx			•	* *		
	Styletto Connect						* *
BTE	Motion C&G X			•	*		

\* Compatible    • Default (Recommended)

Please note – While the X hearing aids will work with the Inductive Charger I, the charger was not custom made for their shape.  
Please verify that the lights turn on when placing the hearing aid in the charger so that you can ensure proper placement.

## Charger Compatibility

Chargers							
	Dry&Clean Charger P X	Dry&Clean Charger SP X	Motion C&G Charger P X	Motion C&G Charger SP X	Insio Charge&Go AX Charger	Active X Pro Charger	Charging Pad
<i>type of charger</i>	Galvanic	Galvanic	Galvanic	Galvanic	MR Inductive	Galvanic	–
<i>part number</i>	10991638	10991644	10996335	10996340	10997503	10991347	purchased separately
							
<i>cord connection</i>	USB-C	USB-C	Attached	Attached	USB-C	USB-C	Attached
Hearing Aids							
BTE	Motion C&G P X	**		*			
	Motion C&G SP X		**		*		
Custom	Insio C&G AX				**		
Misc	Active X Pro					**	*

\* Compatible \*\* Default (Recommended)

Please note – While the X hearing aids will work with the Inductive Charger I, the charger was not custom made for their shape.

Please verify that the lights turn on when placing the hearing aid in the charger so that you can ensure proper placement.

## Rechargeability



### Pure Charge&Go T AX Portable Charger

- Compatible with Pure Charge&Go T AX and CROS Charge&Go AX.
- Portable, closed-lid charger.
- Holds three charges



### Insio Charge&Go AX Charger

- Compatible with Insio C&G AX contactless charging



### Closed-lid Inductive Charger

- Compatible with Pure Charge&Go X and Motion Charge&Go X.
- Dehumidifies, fits custom molds.



### Styletto X Charger

- Compatible with Styletto X.
- Portable, closed-lid charger equipped with Qi Wireless Technology — enables wireless charging on any standard charging pad, including those built into cars and cafe counters.\*



### Dry&Clean Charger

- Dry&Clean Chargers are available for Motion C&G P, Motion C&G SP, and Pure C&G T AX. The world's first hearing aid chargers that simultaneously dry and clean molds/earpieces in just 15 minutes using UV light (UV-C).

\*Wireless charging pad not included with purchase

## Portable AX Charger

For Pure Charge&Go T AX, CROS Charge&Go AX



### NOTICE

Please make sure that the power plug is easily accessible to remove it from power supply if necessary.

- Connect the charger to the power supply. All LEDs for each charging slot turn on successively indicate the charger has power and is ready for use.
- Put the hearing aids in the charging slots as shown in the picture. You can charge both hearing aids at the same time. The hearing aids are automatically turned off and start charging. The charging status is indicated, see table "Charging status".
- When you want to use your hearing aids, take them out of the charger. They will automatically turn on.



## Charger Lights

---

### LED when charger is charging:

---

**Pulsing green:** internal battery is being charged.

---

**Solid green:** charger is fully charged.

---

**Blinking red:** charger error. Please refer to section "Troubleshooting".

---

---

### LED when you are mobile:

---

**Solid green:** the internal battery has enough capacity to fully charge a pair of hearing aids.

---

**Solid orange:** the internal battery does not have enough capacity to fully charge a pair of hearing aids.

---

**Blinking orange:** the internal charger battery is empty. Charging of hearing aids is not possible.

---

## Insio Charge&Go Charger

For Insio Charge&Go AX



### NOTICE

Please make sure that the power plug is easily accessible to remove it from power supply if necessary.

- Connect the charger to the power supply. All LEDs for each charging slot turn on successively indicate the charger has power and is ready for use.
- Put the hearing aids in the charging slots as shown in the picture. You can charge both hearing aids at the same time. The hearing aids are automatically turned off and start charging. The charging status is indicated, see table "Charging status".
- When you want to use your hearing aids, take them out of the charger. They will automatically turn on.



## Charging status

LED	Status description
	<b>No LED:</b> Charger is not connected to power supply.
	Single <b>green</b> LED turns on successively: Charger is connected to a power supply and ready to charge. No hearing aid is inserted in the respective charging slot or the hearing aid is not inserted correctly.
	Flashing <b>orange</b> LED: Charging status of the hearing aid is read out before charging starts.
	Single <b>orange</b> LED: Charging status 0 to 33 %
	Turn on single <b>green</b> LED: Charging status 34 to 66 %
	Turn on two <b>green</b> LEDs: Charging status 67 to 99 %
	Three <b>green</b> LEDs: Hearing aid is fully charged.
	Single <b>orange</b> LED flashes successively: Cooling phase Charging is resumed after 5 minutes.
	Flashing <b>red</b> LEDs: Error Read more in section "Troubleshooting".

## Closed-lid Inductive Charger

For Pure Charge&Go X, CROS Charge&Go X, Motion Charge&Go X BTE



### NOTICE

Please make sure that the power plug is easily accessible to remove it from power supply if necessary.

- Connect the charger to the power supply. All LEDs for each charging slot turn on successively indicate the charger has power and is ready for use.
- Put the hearing aids in the charging slots as shown in the picture. You can charge both hearing aids at the same time. The hearing aids are automatically turned off and start charging. The charging status is indicated, see table "Charging status".
- When you want to use your hearing aids, take them out of the charger. They will automatically turn on.




## Charging status

LED	Status description
	<b>No LED:</b> Charger is not connected to power supply.
	Single <b>green</b> LED turns on successively: Charger is connected to a power supply and ready to charge. No hearing aid is inserted in the respective charging slot or the hearing aid is not inserted correctly.
	Flashing <b>orange</b> LED: Charging status of the hearing aid is read out before charging starts.
	Single <b>orange</b> LED: Charging status 0 to 33 %
	Turn on single <b>green</b> LED: Charging status 34 to 66 %
	Turn on two <b>green</b> LEDs: Charging status 67 to 99 %
	Three <b>green</b> LEDs: Hearing aid is fully charged.
	Single <b>orange</b> LED flashes successively: Cooling phase Charging is resumed after 5 minutes.
	Flashing <b>red</b> LEDs: Error Read more in section "Troubleshooting".

## Stylecto X Charger

- Place the right hearing aid in the charging slot marked with an "R", and the left hearing aid in the charging slot marked with an "L".
- Stow the ear pieces nicely.
- Close the charger's lid for battery power saving reasons.
- The hearing aids are automatically turned off and start charging. The charging status is shown for 7 seconds.



If the status LEDs do not light up when you insert the hearing aids, check the following:

- The hearing aids are not correctly placed. Push them slightly into their charging slots or close the charger's lid.
- The charger is turned off. Press and hold the button for 5 seconds to turn on the charger.
- The charger's internal battery is too low for mobile charging. Connect the charger to an electrical outlet.



**NOTICE:** If you put a hearing aid in the wrong charging slot, the receiver cable may be damaged when closing the charger's lid.

### Charging status of hearing aids

The LED next to the hearing aid indicates the charging status of the hearing aid. Note that the status is only displayed immediately after putting the hearing aids in the charger, when you briefly press the button or when you open the lid. When the charger is off, no status is displayed at all.

#### Hearing Aid LED Status description

**Pulsing green:** Hearing aid is being charged.

**Green:** Hearing aid is fully charged. Internal battery has enough capacity to fully charge two hearing aids.

### Charger LED Status description

**LED is off:** No hearing aid is in the charger, hearing aid is not placed correctly, or charger is off.

**Pulsing orange:** Cooling phase. Charging is resumed when the device has cooled down. Please refer to section "Troubleshooting".

**Blinking red:** Hearing aid error. Please refer to section "Troubleshooting".

**Green:** Internal battery has enough capacity to fully charge two hearing aids.

## **NEW** Dry&Clean Charger

For Pure Charge&Go T AX, Motion P X and Motion SP X



### **Cleaning and Drying**

- Drying and cleaning focuses on tubes, ear tips and molds.
- Put the hearing aids into the charging slots and close the lid.
- UVC cleaning starts automatically after closing the lid and lasts 15 minutes.
- Ongoing UVC cleaning process is indicated by LEDs (see table "Charging and cleaning status").
- Drying starts automatically after closing the lid and lasts 4 hours.


### **Charging**

A full charging cycle takes approximately 3 to 4 hours.





After 30 minutes, it is possible to use your hearing aids for up to 7 hours.




Charging times may vary due to the use of the hearing aids and the age of the power cell.

## Dry&Clean Charger







### CHARGE INDICATORS



Solid yellow 	The charger is detecting the hearing aid.
Flashing green 	The hearing aid is charging, low battery, less than 40%.
Glowing green 	The hearing aid is charging, battery more than 40%.
Solid green 	The hearing aid is fully charged, battery 100% charged.

One indicator flashing red 	An error occurred during charging, and the charging has stopped.
Both indicators flashing red 	An error occurred during charging, and the charging has stopped.
Flashing yellow 	The temperature of the hearing aid has become too high during charging, and charging is paused.

### DRY AND CLEAN INDICATORS



Solid yellow 	The maintenance cycle is paused, but the charging cycle continues.
Flashing white 	The maintenance cycle is in progress (drying/UV-C).
Solid white 	The maintenance cycle is complete.

Flashing red 	An error occurred during the maintenance, the cycle has stopped.
Flashing yellow 	The maintenance cycle is paused due to overheating (fan still running).

## How to order a Charger in ROES

Choose from Other Commodities, and follow these ordering steps:

ROES Main Menu – Patient Orders

**Select an Action:** Other Commodities  
Accessories

**In Commodity Group pull-down menu, choose:** SIVANTOS

**In Vendor pull-down menu, choose:**

In Search box please enter one of the following words:

ROES Name	Description	Part Number
Insio Charge&Go AX Charger	Charger for Insio Charge7Go AX	10997503
Pure Charge&Go AX Portable Charger	Compatible with Pure Charge&Go T AX, CROS Pure AX	10993118
Inductive Charger II	Charger for Pure Charge&Go X, CROS Pure X, and Motion Charge&Go X	10966942
Charger Styletto X	Styletto X charger	10972015
Dry & Clean Charger Motion X BTE P	Dry & Clean Charger for Motion Charge&Go P X BTE	10991638
Dry & Clean Charger Motion X BTE SP	Dry & Clean Charger For Motion Charge&Go SP X BTE	10991644
Smart Lithium-ion Dryer	NX and older charging options	10960141

**NOTE** If you want it delivered to your station, then change the “Deliver To” shipping address, otherwise, it will ship to the patient by default. Please note: a charging unit comes shipped with every Li-ion hearing aid order.

## Portable Power Supply

- The Portable 2600mAh Mobile Power Bank is the perfect power solution when on the road or stuck somewhere without convenient electricity.
- It fits comfortably in a pocket, purse, backpack, or laptop bag to keep Pure Charge&Go or Motion Charge&Go hearing aids charged when needed.
- Capacity of 2600mAh will provide three full charges.






Portable Power  
Supply (10948100)







# RIC Earmolds

Bold text indicates default options

	RIC Earmold 3.0		RIC Embedded Earmold		RIC HP Earmold
					
<b>Styles</b>	Full Shell Half Shell Canal	Canal Lock Skeleton Non-Occluding	Full Shell Half Shell Canal	Canal Lock Skeleton Non-Occluding	Canal
<b>Tubing</b>	Receiver 3.0 S (Size 0-4) Receiver 3.0 M (Size 0-4) Receiver 3.0 P (Size 0-4)	<b>Size 2 is default</b>	Receiver 3.0 S (Size 0-4) Receiver 3.0 M (Size 0-4) Receiver 3.0 P (Size 0-4)	<b>Size 2 is default</b>	Receiver 3.0 HP (Size 0-4) <b>Size 2 is default</b>
<b>Materials</b>	<b>Hard</b>	Soft	<b>Hard</b>	Soft	<b>Hard</b>
<b>Colors</b>	<u>Gloss or Matte Finish</u> Clear, Rose, Beige, Brown	<b>Clear Gloss is default</b>	<u>Gloss or Matte Finish</u> Clear, Rose, Beige, Brown	<b>Clear Gloss is default</b>	Clear, Tan, Medium Brown, Dark Brown
<b>Venting Type and Size</b>	<b>Standard</b> No Vent Semi IROS Trench Vent Open No Vent MicroVent (1.0mm)	Pressure (1.3mm) Small (1.6mm) Medium Select a Vent (2.0mm) Large Select a Vent (3.0mm) X Large (3.5mm) As Large as Possible	<b>Standard</b> No Vent Semi-IROS Trench Vent (Acrylic Only) Pressure Vent (1.0 mm) Small vent (1.6 mm) Medium Vent (2.0 mm)	SAV 2.0 SAV 3.0 SAV 4.0 Large Vent (2.5 mm) XLarge (3.0 mm) As Large as Possible	<b>Standard</b> No Vent  No Vent Micro Vent (1.0 mm) Pressure Vent (1.3mm) <b>Small Vent (1.6mm)</b>
<b>Ear Canal Length</b>	Short, Medium, <b>Long</b>		Short, Medium, <b>Long</b>		Short, Medium, <b>Long</b>
<b>Indicators</b>	<b>Red/Blue Dots</b>	Red/Blue Letters	<b>Red/Blue Dots</b>	Red/Blue Letters	<b>Red/Blue Dots</b> Red/Blue Letters
<b>Dexterity</b>	<b>Removal Line</b>		<b>Removal Line</b>		<b>Removal Line</b>

# BTE Earmolds

Bold text indicates default options

	BTE THIN TUBE 3.0 Earmold		BTE Conventional Earmold	
				
Styles	Full Shell Half Shell Canal	Canal Lock Skeleton Non-Occcluding	Full Shell Half Shell Canal	Canal Lock Skeleton Non-Occcluding
Tubing	Thin Tube 3.0 (Size 0-4) Thin Tube 3.0 Power (Size 0 -4)	<b>Size 2 is default</b>	<b>Standard</b> Thick Thick Dry Extra Thick Libby Horn	Standard Double Bend Clear Thick Double Bend Clear Thick Dry Double Bend Clear Extra Thick Double Bend Clear
Materials	<b>Hard</b>	Soft	<b>Hard</b>	Soft
Colors (See standard color images on next page)	<u>Gloss and Matte Finish Options</u> Clear, Rose, Beige, Brown, Red, Blue	<b>Clear Gloss is default</b>	<u>Gloss and Matte Finish Options</u> Clear, Rose, Beige, Brown, Red, Blue	<b>Clear Gloss is default</b>
Venting Types and Size	Open Pressure Vent (1.0 mm) Small vent (1.6 mm) <b>Medium Vent (2.0 mm)</b>	Large Vent (2.5 mm) XLarge (3.0 mm) As Large as Possible	<b>Standard</b> No Vent  Pressure Vent (1.0 mm) Small vent (1.6 mm) <b>Medium Vent (2.0 mm)</b> SAV 2.0 SAV 3.0	Semi-IROS Trench Vent (Acrylic Only)  SAV 4.0 Large Vent (2.5 mm) XLarge (3.0 mm) As Large as Possible
Ear Canal Length	Short, Medium, <b>Long</b>		Short, Medium, <b>Long</b>	
Indicators	<b>Red/Blue Dots</b>	Red/Blue Letters	<b>Red/Blue Dots</b>	Red/Blue Letters
Dexterity	<b>Removal Line</b>		<b>Removal Line</b>	





## ROES codes

To place an order for any of the following in ROES first choose **Commodities**, then **Parts and Accessories**, then search by item number.

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
<b>RECEIVERS</b>					
MINI RECEIVER 3.0 KIT	MINI RECEIVER 3.0 KIT	10996321	1	\$199.20	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 45 dB SIZE 0 RIGHT EAR	MINIRECEIVER 3.0 45 dB SIZE 0 RIGHT EAR	10992548	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 45 dB SIZE 0 LEFT EAR	MINIRECEIVER 3.0 45 dB SIZE 0 LEFT EAR	10992549	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 45 dB SIZE 1 RIGHT EAR	MINIRECEIVER 3.0 45 dB SIZE 1 RIGHT EAR	10992550	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 45 dB SIZE 1 LEFT EAR	MINIRECEIVER 3.0 45 dB SIZE 1 LEFT EAR	10992552	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 45 dB SIZE 2 RIGHT EAR	MINIRECEIVER 3.0 45 dB SIZE 2 RIGHT EAR	10992553	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 45 dB SIZE 2 LEFT EAR	MINIRECEIVER 3.0 45 dB SIZE 2 LEFT EAR	10992554	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 45 dB SIZE 3 RIGHT EAR	MINIRECEIVER 3.0 45 dB SIZE 3 RIGHT EAR	10992555	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 45 dB SIZE 3 LEFT EAR	MINIRECEIVER 3.0 45 dB SIZE 3 LEFT EAR	10992556	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 45 dB SIZE 4 RIGHT EAR	MINIRECEIVER 3.0 45 dB SIZE 4 RIGHT EAR	10992557	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 45 dB SIZE 4 LEFT EAR	MINIRECEIVER 3.0 45 dB SIZE 4 LEFT EAR	10992558	1	\$51.36	PURE CNG AX, CROS CNG AX

## ROES codes (continued)

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
MINIRECEIVER 3.0 60 dB SIZE 0 RIGHT EAR	MINIRECEIVER 3.0 60 dB SIZE 0 RIGHT EAR	10992559	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 60 dB SIZE 0 LEFT EAR	MINIRECEIVER 3.0 60 dB SIZE 0 LEFT EAR	10992560	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 60 dB SIZE 1 RIGHT EAR	MINIRECEIVER 3.0 60 dB SIZE 1 RIGHT EAR	10992561	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 60 dB SIZE 1 LEFT EAR	MINIRECEIVER 3.0 60 dB SIZE 1 LEFT EAR	10992562	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 60 dB SIZE 2 RIGHT EAR	MINIRECEIVER 3.0 60 dB SIZE 2 RIGHT EAR	10992563	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 60 dB SIZE 2 LEFT EAR	MINIRECEIVER 3.0 60 dB SIZE 2 LEFT EAR	10992564	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 60 dB SIZE 3 RIGHT EAR	MINIRECEIVER 3.0 60 dB SIZE 3 RIGHT EAR	10992565	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 60 dB SIZE 3 LEFT EAR	MINIRECEIVER 3.0 60 dB SIZE 3 LEFT EAR	10992566	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 60 dB SIZE 4 RIGHT EAR	MINIRECEIVER 3.0 60 dB SIZE 4 RIGHT EAR	10992567	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 60 dB SIZE 4 LEFT EAR	MINIRECEIVER 3.0 60 dB SIZE 4 LEFT EAR	10992568	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 65 dB SIZE 0 RIGHT EAR	MINIRECEIVER 3.0 65 dB SIZE 0 RIGHT EAR	10992569	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 65 dB SIZE 0 LEFT EAR	MINIRECEIVER 3.0 65 dB SIZE 0 LEFT EAR	10992570	1	\$51.36	PURE CNG AX, CROS CNG AX

## ROES codes (continued)

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
MINIRECEIVER 3.0 65 dB SIZE 1 RIGHT EAR	MINIRECEIVER 3.0 65 dB SIZE 1 RIGHT EAR	10992571	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 65 dB SIZE 1 LEFT EAR	MINIRECEIVER 3.0 65 dB SIZE 1 LEFT EAR	10992572	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 65 dB SIZE 2 RIGHT EAR	MINIRECEIVER 3.0 65 dB SIZE 2 RIGHT EAR	10992573	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 65 dB SIZE 2 LEFT EAR	MINIRECEIVER 3.0 65 dB SIZE 2 LEFT EAR	10992574	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 65 dB SIZE 3 RIGHT EAR	MINIRECEIVER 3.0 65 dB SIZE 3 RIGHT EAR	10992575	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 65 dB SIZE 3 LEFT EAR	MINIRECEIVER 3.0 65 dB SIZE 3 LEFT EAR	10992576	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 65 dB SIZE 4 RIGHT EAR	MINIRECEIVER 3.0 65 dB SIZE 4 RIGHT EAR	10992577	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER 3.0 65 dB SIZE 4 LEFT EAR	MINIRECEIVER 3.0 65 dB SIZE 4 LEFT EAR	10992578	1	\$51.36	PURE CNG AX, CROS CNG AX
MINIRECEIVER COMPLETE SET	MINI RECEIVER 2.0 KIT	10960146	1	\$546.00	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
S-RECEIVER 2.0 SIZE 0 L	MINI RECEIVER 2.0 45 dB SIZE 0 LEFT EAR	10669907	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
S-RECEIVER 2.0 SIZE 0 R	MINI RECEIVER 2.0 45 dB SIZE 0 RIGHT EAR	10669908	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
S-RECEIVER 2.0 SIZE 1 L	MINI RECEIVER 2.0 45 dB SIZE 1 LEFT EAR	10669909	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC

## ROES codes (continued)

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
S-RECEIVER 2.0 SIZE 1 R	MINI RECEIVER 2.0 45 dB SIZE 1 RIGHT EAR	10669910	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
S-RECEIVER 2.0 SIZE 2 L	MINI RECEIVER 2.0 45 dB SIZE 2 LEFT EAR	10669911	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
S-RECEIVER 2.0 SIZE 2 R	MINI RECEIVER 2.0 45 dB SIZE 2 RIGHT EAR	10669912	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
S-RECEIVER 2.0 SIZE 3 L	MINI RECEIVER 2.0 45 dB SIZE 3 LEFT EAR	10669913	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
S-RECEIVER 2.0 SIZE 3 R	MINI RECEIVER 2.0 45 dB SIZE 3 RIGHT EAR	10669914	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
M-RECEIVER 2.0 SIZE 0 L	MINI RECEIVER 2.0 60 dB SIZE 0 LEFT EAR	10824820	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
M-RECEIVER 2.0 SIZE 0 R	MINI RECEIVER 2.0 60 dB SIZE 0 RIGHT EAR	10824821	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
M-RECEIVER 2.0 SIZE 1 L	MINI RECEIVER 2.0 60 dB SIZE 1 LEFT EAR	10824822	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
M-RECEIVER 2.0 SIZE 1 R	MINI RECEIVER 2.0 60 dB SIZE 1 RIGHT EAR	10824823	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
M-RECEIVER 2.0 SIZE 2 L	MINI RECEIVER 2.0 60 dB SIZE 2 LEFT EAR	10824824	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
M-RECEIVER 2.0 SIZE 2 R	MINI RECEIVER 2.0 60 dB SIZE 2 RIGHT EAR	10824825	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
M-RECEIVER 2.0 SIZE 3 L	MINI RECEIVER 2.0 60 dB SIZE 3 LEFT EAR	10824826	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC

## ROES codes (continued)

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
M-RECEIVER 2.0 SIZE 3 R	MINI RECEIVER 2.0 60dB SIZE 3 RIGHT EAR	10824827	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
M-RECEIVER 2.0 SIZE 4 L	MINI RECEIVER 2.0 60 dB SIZE 1 LEFT EAR	10933054	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
M-RECEIVER 2.0 SIZE 4 R	MINI RECEIVER 2.0 60 dB SIZE 1 RIGHT EAR	10933053	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
P-RECEIVER 2.0 SIZE 0 L	MINI RECEIVER 2.0 70 dB SIZE 0 LEFT EAR	10669923	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
P-RECEIVER 2.0 SIZE 0 R	MINI RECEIVER 2.0 70 dB SIZE 0 RIGHT EAR	10669924	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
P-RECEIVER 2.0 SIZE 1 L	MINI RECEIVER 2.0 70 dB SIZE 1 LEFT EAR	10669925	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
P-RECEIVER 2.0 SIZE 1 R	MINI RECEIVER 2.0 70 dB SIZE 1 RIGHT EAR	10669926	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
P-RECEIVER 2.0 SIZE 2 L	MINI RECEIVER 2.0 70 dB SIZE 2 LEFT EAR	10669927	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
P-RECEIVER 2.0 SIZE 2 R	MINI RECEIVER 2.0 70 dB SIZE 2 RIGHT EAR	10669928	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
P-RECEIVER 2.0 SIZE 3 L	MINI RECEIVER 2.0 70 dB SIZE 3 LEFT EAR	10669929	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
P-RECEIVER 2.0 SIZE 3 R	MINI RECEIVER 2.0 70 dB SIZE 3 RIGHT EAR	10669930	1	\$51.36	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
<b>THIN TUBES</b>					
3.0 THIN TUBE KIT	THIN TUBE 3.0 FITTING KIT	10992763	1	\$93.60	MOTION CNG X BTE, MOTION CNG P X BTE

## ROES codes (continued)

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
BINAURAL THINTUBE 3.0 S0	THIN TUBE 3.0 SIZE 0	10991609	2	\$8.53	MOTION CNG X BTE , MOTION CNG P X BTE
BINAURAL THINTUBE 3.0 S1	THIN TUBE 3.0 SIZE 1	10991610	2	\$8.53	MOTION CNG X BTE , MOTION CNG P X BTE
BINAURAL THINTUBE 3.0 S2	THIN TUBE 3.0 SIZE 2	10991611	2	\$8.53	MOTION CNG X BTE , MOTION CNG P X BTE
BINAURAL THINTUBE 3.0 S3	THIN TUBE 3.0 SIZE 3	10991612	2	\$8.53	MOTION CNG X BTE , MOTION CNG P X BTE
BINAURAL THINTUBE 3.0 S4	THIN TUBE 3.0 SIZE 4	10991613	2	\$8.53	MOTION CNG X BTE , MOTION CNG P X BTE
BINAURAL THINTUBE 3.0 S0 POWER	POWER THIN TUBE 3.0 SIZE 0	10992492	2	\$8.53	MOTION CNG X BTE , MOTION CNG P X BTE
BINAURAL THINTUBE 3.0 S1 POWER	POWER THIN TUBE 3.0 SIZE 1	10992493	2	\$8.53	MOTION CNG X BTE , MOTION CNG P X BTE
BINAURAL THINTUBE 3.0 S2 POWER	POWER THIN TUBE 3.0 SIZE 2	10992494	2	\$8.53	MOTION CNG X BTE , MOTION CNG P X BTE
BINAURAL THINTUBE 3.0 S3 POWER	POWER THIN TUBE 3.0 SIZE 3	10992495	2	\$8.53	MOTION CNG X BTE , MOTION CNG P X BTE
BINAURAL THINTUBE 3.0 S4 POWER	POWER THIN TUBE 3.0 SIZE 4	10992496	2	\$8.53	MOTION CNG X BTE , MOTION CNG P X BTE
TUBE THIN 2.0 FITTING SET	THIN TUBE 2.0 FITTING KIT	10828299	1	\$169.71	MOTION CNG NX BTE, MOTION 13 NX BTE, MOTION 13P NX BTE
THIN TUBE 2.0 1L	THIN TUBE 2.0 SIZE 1 LEFT	10828291	5	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE, MOTION 13P NX BTE
THIN TUBE 2.0 1R	THIN TUBE 2.0 SIZE 1 RIGHT	10828292	5	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE, MOTION 13P NX BTE
THIN TUBE 2.0 2L	THIN TUBE 2.0 SIZE 2 LEFT	10828293	5	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE, MOTION 13P NX BTE

## ROES codes (continued)

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
THIN TUBE 2.0 2R	THIN TUBE 2.0 SIZE 2 RIGHT	10828294	5	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE, MOTION 13P NX BTE
THIN TUBE 2.0 3L	THIN TUBE 2.0 SIZE 3 LEFT	10828295	5	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE, MOTION 13P NX BTE
THIN TUBE 2.0 3R	THIN TUBE 2.0 SIZE 3 RIGHT	10828296	5	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE, MOTION 13P NX BTE
THIN TUBE 2.0 4L	THIN TUBE 2.0 SIZE 4 LEFT	10828297	5	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE, MOTION 13P NX BTE
THIN TUBE 2.0 4R	THIN TUBE 2.0 SIZE 4 RIGHT	10828298	5	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE, MOTION 13P NX BTE
<b>DOMES AND SLEEVES</b>					
EARTIP 3.0 5MM OPEN	EARTIP 3.0 S OPEN	10993678	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
EARTIP 3.0 7MM OPEN	EARTIP 3.0 M OPEN	10993680	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
EARTIP 3.0 10MM OPEN	EARTIP 3.0 L OPEN	10993681	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
EARTIP 3.0 8MM TULIP	EARTIP 3.0 M TULIP	10993682	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
EARTIP 3.0 12MM TULIP	EARTIP 3.0 L TULIP	10993683	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
SLEEVE 3.0 VENTED XS	SLEEVE 3.0 XS VENTED	10992863	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
SLEEVE 3.0 VENTED S	SLEEVE 3.0 VENTED S TBLK	10992864	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE

## ROES codes (continued)

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
SLEEVE 3.0 VENTED M	SLEEVE 3.0 VENTED M TBLK	10992865	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
SLEEVE 3.0 VENTED L	SLEEVE 3.0 VENTED L TBLK	10992866	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
SLEEVE 3.0 CLOSED XS	SLEEVE 3.0 CLOSED XS TBLK	10992867	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
SLEEVE 3.0 CLOSED S	SLEEVE 3.0 CLOSED S TBLK	10992868	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
SLEEVE 3.0 CLOSED M	SLEEVE 3.0 CLOSED M TBLK	10992869	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
SLEEVE 3.0 CLOSED L	SLEEVE 3.0 CLOSED L TBLK	10992870	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
SLEEVE 3.0 POWER XS	SLEEVE 3.0 POWER XS CXX	10997192	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
SLEEVE 3.0 POWER S	SLEEVE 3.0 POWER S CXX	10997193	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
SLEEVE 3.0 POWER M	SLEEVE 3.0 POWER M CXX	10997194	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
SLEEVE 3.0 POWER L	SLEEVE 3.0 POWER L CXX	10997195	6	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P BTE
CLICK SLEEVES XS VENTED	CLICK SLEEVE 2.0 XS VENTED	10951152	6	\$8.53	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
CLICK SLEEVES S VENTED	CLICK SLEEVE 2.0 S VENTED	10951146	6	\$8.53	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC

## ROES codes (continued)

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
CLICK SLEEVES M VENTED	CLICK SLEEVE 2.0 M VENTED	10951148	6	\$8.53	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
CLICK SLEEVES L VENTED	CLICK SLEEVE 2.0 L VENTED	10951150	6	\$8.53	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
CLICK SLEEVES S CLOSED	CLICK SLEEVE 2.0 S CLOSED	10951147	6	\$8.53	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
CLICK SLEEVES M CLOSED	CLICK SLEEVE 2.0 M CLOSED	10951149	6	\$8.53	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
CLICK SLEEVES L CLOSED	CLICK SLEEVE 2.0 L CLOSED	10951151	6	\$8.53	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
CLICK DOMES OPEN 4MM	CLICK DOME S OPEN	10426016	6	\$5.12	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
CLICK DOMES OPEN 8MM	CLICK DOME M OPEN	10426014	6	\$5.12	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
CLICK DOMES OPEN 10MM	CLICK DOME L OPEN	10426015	6	\$5.12	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
CLICK DOMES CLOSED 6MM	CLICK DOME S CLOSED	10426017	6	\$5.12	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
CLICK DOMES CLOSED 8MM	CLICK DOME M CLOSED	10426012	6	\$5.12	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
CLICK DOMES CLOSED 10MM	CLICK DOME L CLOSED	10426013	6	\$5.12	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
CLICK DOMES DOUBLE 8/10MM	CLICK DOME S DOUBLE	10426018	6	\$5.12	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC

## ROES codes (continued)

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
CLICK DOMES DOUBLE 10/12MM	CLICK DOME M DOUBLE	10426019	6	\$5.12	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
DOMES 8/12MM SEMI OPEN CLICK RIC	CLICK DOME M SEMI OPEN	10670790	6	\$5.12	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
LIFE TIP 8 MM OPEN	LIFE TIP REFILL 8MM	10054916	1	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE. MOTION 13P NX BTE
LIFE TIP 8 MM CLOSED	LIFE TIP REFILL 8MM FULL	10169330	1	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE. MOTION 13P NX BTE
LIFE TIP 4 MM	LIFE TIP REFILL 4 MM	10169328	1	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE. MOTION 13P NX BTE
LIFE TIP 6 MM	LIFE TIP REFILL 6MM	10169329	1	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE. MOTION 13P NX BTE
LIFE TIP 10 MM	LIFE TIP REFILL 10MM	10054917	1	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE. MOTION 13P NX BTE
DOUBLE TIP 8/10	LIFE TIP DOUBLE REFILL 8/10MM	10175566	1	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE. MOTION 13P NX BTE
DOUBLE TIP 10/12	LIFE TIP DOUBLE REFILL 10/12MM	10175715	1	\$8.53	MOTION CNG NX BTE, MOTION 13 NX BTE. MOTION 13P NX BTE
<b>CHARGERS</b>					
INSIO CNG AX CHARGER	CHARGER FOR INSIO CUSTOM HEARING AIDS	10997503	1	\$175.76	INSIO CNG AX
PORTABLE CHARGER FOR PURE CNG AX	PORTABLE CHARGER FOR PURE C&G AX	10993118	1	\$175.76	PURE CNG AX, CROS CNG AX

## ROES codes (continued)

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
DRY & CLEAN CHARGER FOR PURE CNG AX	DRY&CLEAN CHARGER FOR PURE C&G AX	10991650	1	\$218.40	PURE CNG AX, CROS CNG AX
STANDARD CHARGER FOR PURE CNG AX	CHARGER FOR PURE C&G AX	10996600	1	\$175.76	PURE CNG AX, CROS CNG AX
DRY & CLEAN CHARGER MOTION X BTE P	DRY & CLEAN CHARGER FOR MOTION CNG P X BTE	10991638	1	\$218.40	MOTION CNG X P BTE
DRY & CLEAN CHARGER MOTION X BTE SP	DRY & CLEAN CHARGER FOR MOTION CNG SP X BTE	10991644	1	\$218.40	MOTION CNG X SP BTE
CHARGER MOTION XBTE P	CHARGER FOR MOTION CNG P X BTE	10996335	1	\$175.76	MOTION CNG P BTE
CHARGER MOTION X BTE SP	CHARGER FOR MOTION CNG SP X BTE	10996340	1	\$175.76	MOTION CNG SP BTE
INDUCTIVE CHARGER II	INDUCTIVE CHARGER II FOR PURE AND MOTION	10966942	1	\$165.57	PURE CNG X RIC, MOTION CNG X, MOTION CNG NX
CHARGER STYLETTO X	CHARGER FOR STYLETTO X	10972015	1	\$175.76	STYLETTO X RIC
CHARGER ACTIVE PRO	CHARGER FOR ACTIVE PRO	10991347	1	\$175.76	ACTIVE PRO
CHARGER INDUCTIVE SET LI-ION	LI-ION INDUCTIVE CHARGER SET	10948199	1	\$175.76	CELLION PX, PURE CNG NX, MOTION CNG NX
<b>TELECOIL DOORS AND TOOLS</b>					
TELECOIL DOOR SET BG	TELECOIL DOOR SET BEIGE	10950811	1	\$16.64	Pure 13 T NX RIC, Motion 13 NX BTE, Motion 13 P NX BTE
TELECOIL DOOR SET DGT	TELECOIL DOOR SET DGT	10950812	1	\$16.64	Pure 13 NX RIC, Motion 13 NX BTE, Motion 13 P NX BTE
TELECOIL DOOR SET BLK	TELECOIL DOOR SET BLACK	10953591	1	\$16.64	Pure 13 NX RIC, Motion 13 NX BTE, Motion 13 P NX BTE
TELECOIL DOOR SET DBR	TELECOIL DOOR SET DARK BROWN	10953592	1	\$16.64	Pure 13 NX RIC, Motion 13 NX BTE, Motion 13 P NX BTE

## ROES codes (continued)

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
TELECOIL DOOR SET DKC	TELECOIL DOOR SET DKC	10953764	1	\$16.64	Pure 13 NX RIC, Motion 13 NX BTE, Motion 13 P NX BTE
TELECOIL DOOR SET GB	TELECOIL DOOR SET GB	10953765	1	\$16.64	Pure 13 NX RIC, Motion 13 NX BTE, Motion 13 P NX BTE
TELECOIL DOOR SET GNT	TELECOIL DOOR SET GNT	10953593	1	\$16.64	Pure 13 NX RIC, Motion 13 NX BTE, Motion 13 P NX BTE
TELECOIL DOOR SET GR	TELECOIL DOOR SET GR	10953594	1	\$16.64	Pure 13 NX RIC, Motion 13 NX BTE, Motion 13 P NX BTE
TELECOIL DOOR SET PRL	TELECOIL DOOR SET PRL	10953766	1	\$16.64	Pure 13 NX RIC, Motion 13 NX BTE, Motion 13 P NX BTE
TELECOIL DOOR SET SB	TELECOIL DOOR SET SB	10953767	1	\$16.64	Pure 13 NX RIC, Motion 13 NX BTE, Motion 13 P NX BTE
TELECOIL DOOR SET SLV	TELECOIL DOOR SET SILVER	10953768	1	\$16.64	Pure 13 NX RIC, Motion 13 NX BTE, Motion 13 P NX BTE
MAGNET AUTOPHONE	TWINPHONE MAGNET	10940330	1	\$16.64	Pure 13 NX RIC, Motion 13 NX BTE, Motion 13 P NX BTE
<b>WAXGUARDS</b>					
WAXGUARD 3.0 NANOCARE	EARWEAR 3.0 NANOCARE WAXGUARD	10993607	1	\$0.57	PURE CNG AX, INSIO AX, MOTION CNG X BTE, MOTION CNG P X BTE
WAXGUARD RIC MINI	MINI RECEIVER WAXGUARD	10966257	10	\$2.98	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
QUICKGUARD PACK	CXX QuickGuard	10962504	1	\$3.41	INSIO NX CUSTOMS
WAXGUARD C-GRID	C-GRID WAXGUARD	10423507	16	\$8.53	INSIO NX CUSTOMS
WAX SPRING,RCVR TUBE ITE	WAX SPRING FOR RCVR TUBE	7080208	1	\$0.00	INSIO NX CUSTOMS
<b>CONCHA LOCKS</b>					
CONCHA LOCK 3.0	EARWEAR 3.0 CONCHA LOCK	10993677	10	\$8.53	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P X BTE
RECEIVER CONCHA LOCK S45dB	MINI RECEIVER 2.0 CONCHA LOCK FOR S RECEIVER	10601877	10	\$0.00	PURE CNG X RIC, PURE 312T X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC

## ROES codes (continued)

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
RECEIVER CONCHA LOCK M60dB	MINI RECEIVER 2.0 CONCHA LOCK FOR M RECEIVER	10937826	10	\$0.00	PURE CNG X RIC, PURE 312T X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
RECEIVER CONCHA LOCK P70dB	MINI RECEIVER 2.0 CONCHA LOCK FOR P RECEIVER	10673516	10	\$0.00	PURE CNG X RIC, PURE 312T X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
<b>EARHOOKS</b>					
EARHOOK MOTION C-G NX	EARHOOK FOR MOTION CNG NX BTE	10959864	1	\$0.00	MOTION CNG NX BTE
MOTION PX/NX EARHOOK	EARHOOK FOR MOTION 13 NX BTE	10947168	1	\$0.00	MOTION 13 NX BTE
THINTUBE ADAPTER MOTION	THIN TUBE ADAPTER FOR MOTION CNG NX BTE	10958343	1	\$3.12	MOTION CNG NX BTE
MOTION PX/NX TUBE ADAPTER	THIN TUBE ADAPTER FOR MOTION 13 NX BTE	10947303	1	\$0.00	MOTION 13 NX BTE
EARHOOK MOTION X SMALL	EARHOOK MOTION X SMALL FOR	10980403	1	\$0.00	MOTION CNG X BTE, MOTION CNG P X BTE
EARHOOK MOTION X STANDARD	EARHOOK MOTION X STANDARD	10980402	1	\$0.00	MOTION CNG X BTE, MOTION CNG P X BTE
EARHOOK MOTION X SP SMALL	EARHOOK MOTION X SP SMALL	10993620	1	\$0.00	MOTION CNG SP X BTE
<b>MIC COVERS</b>					
CUSTOM MIC COVER BEIGE	MIC COVERS FOR CUSTOM HEARING AIDS - BEIGE	10541504	10	\$0.00	INSIO C&G AX, INSIO NX
CUSTOM MIC COVER BROWN	MIC COVER FOR CUSTOM HEARING AIDS - BROWN	10541424	1	\$0.00	INSIO C&G AX, INSIO NX
CUSTOM MIC COVER MOCHA	MIC COVER FOR CUSTOM HEARING AIDS - MOCHA	10541423	1	\$0.00	INSIO C&G AX, INSIO NX

## ROES codes (continued)

VA ROES ITEM	DESCRIPTION	ITEM #	QUANTITY	PRICE	COMPATIBLE PRODUCTS
CUSTOM MIC COVER DARK BROWN AND BLACK	MIC COVER FOR CUSTOM HEARING AIDS - DARK BROWN AND BLACK	10541574	1	\$0.00	INSIO C&G AX, INSIO NX
<b>MISC TOOLS</b>					
3.0 MEASURING TOOL	EARWEAR 3.0 MEASURING TOOL	10991403	1	\$0.00	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P X BTE
EAR-TIP 3.0 INSERTION TOOL	EARWEAR 3.0 INSERTION TOOL	10994214	1	\$0.00	PURE CNG AX, CROS CNG AX, MOTION CNG X BTE, MOTION CNG P X BTE
TOOL MEASURING	RECEIVER MEASURING TOOL	10174571	1	\$0.00	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
MINI RECEIVER TOOL	MINI RECEIVER REMOVAL TOOL	10674370	1	\$0.00	PURE CNG X, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX, PURE 312T NX RIC, PURE 13 NX RIC
THIN TUBE 3.0 PIN EJECTOR TOOL	EARHOOK REMOVAL TOOL	10600047	1	\$0.00	MOTION CNG X BTE, MOTION CNG P X BTE
TOOL CLICK MOLD	CLICK MOLD SERVICE TOOL	10597930	1	\$0.00	PURE CNG X RIC, PURE 312T X RIC, STYLETTO X RIC, PURE CNG NX RIC, PURE 312T NX RIC, PURE 13 NX RIC
EARTIP 3.0 HANDLING TOOL	3.0 EARMOLD SERVICE TOOL	10984033	1	\$0.00	MOTION CNG X BTE, MOTION CNG P X BTE
MASK EXTENDERS	MASK EXTENDERS	10242249	1	\$0.00	ALL MODELS
<b>ACCESSORY PARTS</b>					
STREAMLINE TV POWER SUPPLY	POWER SUPPLY FOR STREAMLINE TV	10944815	1	\$0.00	STREAMLINE TV
SMARTCONNECT POWER ADAPT SET	SMARTCONNECT POWER ADAPT SET	10824346	1	\$0.00	STREAMLINE TV

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners.

Android is a trademark of Google Inc.

iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.

"Made for iPhone" means that an electronic accessory has been designed to connect specifically to iPhone and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPhone may affect wireless performance.



**Our mission is to help improve the quality of life of Veterans,  
Active Military and American Indians with better hearing.**

**Signia Government Services: (800) 955-4327 | Consumer Hotline: (800) 350-6093**