

Powered by the Signia Augmented Xperience Platform



Expanding on our mission of enhancing human performance

Signia's latest innovation – the Pure Charge&Go AX – is powered by the NEW, world's first Augmented Xperience platform which delivers up to 96 channels with two dedicated processors and a host of clinically proven benefits including NEW and Enhanced audiological features. Ultimately the NEW Augmented Xperience platform is designed to deliver outstanding speech quality in quiet, moderate, and noisy environments!

Provider Benefits

- Fully featured RIC designed to **support a broad range of patient's needs:**
 - Easily connect to iPhone and ASHA compatible Android devices**
 - Integrated **T-coil in a compact RIC design and CROS** options
 - Rocker switch for onboard control
 - Sleek design combined with new EarWear 3.0 which is proven to deliver higher physical fit rates and better retention for increased patient satisfaction
- Rechargeable with **multiple charger options** to choose from to appeal to each patient's unique needs
 - Portable Charger** – provides **3 days of extra charges** for patients on the go
 - Pure Charger** – simplifies charging with its **small, freestanding design**
 - Dry&Clean Charger** – to **keep hearing aids in top working conditions**

Pure Charger and Dry&Clean Charger are available separately

Wearer Benefits

- Sleek hearing aid design** that is both comfortable and discreet (via EarWear 3.0)
- Easily connects to iPhone and ASHA compatible Android devices** for effortless streaming from any Bluetooth® device.
- Fully featured hearing aid with optional support** of induction looping systems and a CROS transmitter for single-sided deafness
- Multiple charger options to match your lifestyle with **up to 36 hours of wear time on a single charge.**

Key Features*

- Up to 36 hours of use on a single charge with T-coil, and up to 24 hours without T-coil (reflects 5 hours of streaming in both options).
- Powered by the **World's First Augmented Xperience platform with two dedicated processors**
 - More than **25% better speech understanding in noise** compared to own hearing aids¹
 - 95% of participants reported exceptional speech understanding in their home environment²
- Remote Care Ready

* Refer to Feature Overview Chart in the AX Platform Brochure

¹ Signia White Paper. Based on transformation of the 1.8 dB SNR improvement using a value of the slope of the psychometric function (Wagener & Brand, 2005)

² Signia White Paper. 18/19 participants had a mean rating of speech intelligibility-related questions of more than 8.0 and 9.0, respectively



24/7 Digital Assistant



Bluetooth Connectivity



Li-ion Rechargeability



Own Voice Processing



Dynamic Soundscape Processing 2.0



Acoustic Motion Sensors



Signia App



Tinnitus



T-coil

Available Colors:



Signia App is equipped with the Signia Assistant for 24/7 support plus Signia Face Mask Mode Improves audibility for those communicating with someone wearing a face mask - all with the touch of a button in the Signia App.

Type	Pure Charge&Go AX
Available Performance Level	7
Matrix	S (110/46) M (119/60) P (122/65) HP (131/75)
Battery Size	Lithium-ion power cell
Features	
Frequency channels/ compression controls	96/20 (48 channels per processor)
Hearing Programs	6
Bluetooth® Streaming	•
TwinPhone ¹	•
TeleCare	•
Tinnitus therapies	•
CROS/BiCROS ²	•
Rocker Switch ³	•
Rechargeable	•
T-coil	◦
App/Accessory Compatibility	
Signia App	•
Streamline TV/Streamline Mic	•
miniPocket	•
Noahlink Wireless	•
Fitting set	EarWear 3.0

Pure Charge&Go AX works exclusively with the new EarWear 3.0.

Multiple charger options available



Pure Portable Charger
Standard Option for Level
5 AX and 7 AX Hearing Aids



Pure Charger
Standard Option for Level
3 AX Hearing Aids



Pure Dry&Clean Charger
Upgrade Option for Level 3 AX,
5 AX and 7 AX Hearing Aids

Pure Charger and Dry&Clean Charger are available separately

Compatible with:



Streamline™ TV



Streamline™ Mic



miniPocket™

¹ Activate via rocker switch, app or accessory

² Compatible with CROS Pure Charge&Go AX

³ Programmable via Connexx

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners.

- Standard
- Optional