

# Signia Xperience CROS solutions:

Three revolutionary world's first innovations.



The ultimate rechargeable CROS/BiCROS experience

## RECHARGEABLE

### CROS Pure® Charge&Go X Transmitter

- World's first CROS/BiCROS technology to benefit from acoustic-motion sensors for the clearest speech understanding that adapts to the wearer's movements
- Own Voice Processing (OVP™) with BiCROS
- Li-ion rechargeability and auto on/auto off function when removed from charger for easy handling
- Directional pattern skewed toward the unaidable ear for improved speech understanding
- CROS mixing-level adjustment available in the Signia App to balance spatial benefit vs. noise level
- Compatible with the Signia App and miniPocket™



The easy-to-use, direct streaming CROS/BiCROS experience

## SIMPLE & CONNECTED

### CROS Pure® 312 X Transmitter

- World's first CROS/BiCROS technology to benefit from acoustic-motion sensors for the clearest speech understanding that adapts to the wearer's movements
- Own Voice Processing (OVP) with BiCROS
- Directional pattern skewed toward the unaidable ear for improved speech understanding
- CROS mixing level adjustment available with the Signia App to balance spatial benefit vs. noise level
- Compatible with the Signia App and miniPocket

**YOURSOUND** TECHNOLOGY allows you to **HEAR WHAT MATTERS** to you.

**ACOUSTIC-MOTION  
SENSORS**

**DYNAMIC SOUNDSCAPE  
PROCESSING**

**OWN VOICE  
PROCESSING**

# Signia X CROS compatibility chart.



		CROS Pure Charge&Go X	CROS Pure 312 X
Products	<b>Pure Charge&amp;Go X</b> Max Output/Gain: 130/75 dB	★	
	<b>Pure 312 X (T-coil optional)*</b> Max Output/Gain: 130/75 dB		★
Apps & Accessories	<b>StreamLine™ TV</b>	HI	HI
	<b>StreamLine™ Mic</b>	HI	HI
	<b>Signia App</b>	HI + CROS	HI + CROS
	<b>miniPocket</b>	HI	HI

★ Recommended compatibility for ease of patient use

HI - Hearing instrument

\* CROS transmitter mics automatically turn off when telecoil program is active

## Signia App



- Displays features dependent on hearing aid model
- Remote control, streaming, and TeleCare options

## StreamLine TV



- Dolby Digital® stereo TV streaming to hearing aids
- Adjust TV-streamed volume via the Signia App

## StreamLine Mic



- Hands-free streaming of phone calls and music
- Remote mic for unmatched speech clarity in noise

## miniPocket



- Make volume and program changes
- Small enough to carry on a key ring

## CROS Pure Charge&Go X in ear



## CROS Pure 312 X in ear



Dolby and the double-D symbol are registered trademarks of Dolby Laboratories. All other trademarks remain the property of their respective owners.

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice. The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners. Android is a trademark of Google Inc.

Copyright © 2020 Signia GmbH. All rights reserved. 4/20 1.0 SIV-7568-20