

Stylecto X.

The Complete Hearing Experience

Performance levels

7X



reddot design award
winner 2019



NEW
Portable Charger

HIGHLIGHTS:

- Most natural own voice* and best speech understanding in noise**
- High-quality audio streaming
- Ultra HD e2e and Bluetooth connectivity
- Automatic on/off functionality with charger
- Lasts all day: Up to 19 hours of listening time (16 hours with 5 hours of streaming)
- 30-minute fast charge = 5 more hours of use

A fully personalized Xperience providing improved hearing in noise while adapting to your movements and maintaining awareness of the environment around you.



Key features & benefits

- Highly sophisticated hearwear™ in a unique, patient-preferred† design
- Portable charger holds **3 days of charges** and features **integrated Qi wireless technology** for simplified charging
- **Enhanced functionality** to support a wider fitting range with exchangeable receivers



Signia App ready

- **NEW** Signia Assistant
- TeleCare™
- Remote control
- Connectivity & streaming



Equipped with *Bluetooth** for easy and carefree enjoyment of TV, music, and phone calls



Rechargeable



NEW Charger equipped with Qi Wireless Technology - enables wireless charging on any standard charging pad, including those built into cars and cafe counters ††

Available colors



Black / Silver Cosmic Blue/ Rose Gold Snow White/ Rose Gold Black White



StreamLine TV
TV signal is directly streamed into the hearing aids



StreamLine Mic
Hands-free audio streaming (phone calls, music, etc.) from any device with Bluetooth connectivity



miniPocket
Discreet remote control without the app

YOURSOUND TECHNOLOGY allows you to **HEAR WHAT MATTERS** to you.

ACOUSTIC-MOTION SENSORS

DYNAMIC SOUNDSCAPE PROCESSING

OWN VOICE PROCESSING

Type

Styletto

Available performance levels

7X

Picture



Product details

Matrix (2 ccm)

S 45/108 dB
M 60/119 dB
P 70/124 dB

Battery size

Lithium-ion power cell

Power consumption

1.9 - 2.1 mA

Features

Frequency channels/ Compression controls

48/20

Hearing programs

6

Ultra HD e2e wireless

●¹

Direct streaming

●

TwinPhone

●⁵

TeleCare

●

Tinnitus therapies

●

CROS/BiCROS

-

Push button

-

Volume control

●³

Rechargeable

●

T-coil

-

Battery door on/off function

-

Audio input

●⁷

Modifications

Color conversion kit

-

App/Accessory compatibility

Signia App

●

StreamLine TV/StreamLine Mic

●

miniPocket

●

ConnexxAir

-

Noahlink Wireless

●

Programming adapter/cable

-


Fitting set

-

- Standard
- Optional


- (1) Ultra HD e2e with binaural audio signal processing
- (3) To be activated via Signia App or miniPocket
- (5) Requires bilateral fitting
- (7) Via Bluetooth Low Energy

S-Receiver




Open Click Domes
Double Click Domes
Click Mold (no vent)

M-Receiver



Open Click Domes
Double Click Domes
Click Mold (no vent)

P-Receiver



Double Click Domes
Click Mold (no vent)



† Hakvoort, C. & Burton, B. (2019). Increasing Style, Reducing Stigma: The Styletto Solution.

** Wireless charging pad not included with purchase

* 2017 "OVP Study" conducted at University of Northern Colorado examining the effect of Own Voice Processing on spontaneous acceptance among dissatisfied wearers after first fit of hearing aids. Further details: www.signia-pro.com/ovp-study

** Froehlich M, Freels K, Powers T. (2015). Speech recognition benefit obtained from binaural beamforming hearing aids: comparison to omnidirectional and individuals with normal hearing. *AudiologyOnline*, Article 14338. Retrieved from <http://www.audiologyonline.com>

The *Bluetooth*® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners.

The 'Qi' symbol is a trademark of the Wireless Power Consortium.

Copyright © 2020 Signia GmbH. All rights reserved. SH202070 3/20 1.0 V11011980 SIV-8824-20