

Signia Pure Charge&Go T AX

Powered by the Signia Augmented Xperience Platform



Expanding on our mission of enhancing human performance

Signia's Pure Charge&Go T AX – is powered by the world's first Augmented Xperience platform which delivers up to 96 channels with two dedicated processors for split processing and a host of clinically proven benefits including NEW and enhanced audiological features. Ultimately the Augmented Xperience platform is designed to deliver outstanding speech quality in quiet, moderate, and noisy environments!

Provider Benefits

- Fully featured RIC designed to **support a broad range of patient's needs:**
 - Easily connect to iPhone® and ASHA compatible Android™ devices**
 - Integrated **T-coil in a compact RIC design and CROS** options
 - Rocker switch for onboard control
 - Sleek design combined with new EarWear 3.0 which is proven to deliver higher physical fit rates and better retention for increased patient satisfaction
- Rechargeable with **multiple charger options** to choose from to appeal to each patient's unique needs
 - Pure Charger** – simplifies charging with its **small, freestanding design**
 - Portable Charger** – provides **3 days of extra charges** for patients on the go
 - Dry&Clean Charger** – to **keep hearing aids in top working condition**

Wearer Benefits

- Sleek hearing aid design** that is both comfortable and discreet (via EarWear 3.0)
- Easily connects to iPhone and ASHA compatible Android devices** for effortless streaming from any Bluetooth® device.

- Fully featured hearing aid with optional support** of induction looping systems and a CROS transmitter for single-sided deafness
- Multiple charger options to match your lifestyle with **up to 36 hours of wear time on a single charge** with T-coil, and up to 24 hours without T-coil (reflects 5 hours of streaming in both options).
- Remote Care Ready

Key Features*

- Powered by the **World's First Augmented Xperience platform with two dedicated processors**
 - More than **25% better speech understanding in noise** compared to own hearing aids
 - 94% of wearers report excellent speech understanding** with Signia AX*

New Features

- Own Voice Processing (OVP™) 2.0** utilizes split processing of the Augmented Xperience platform to separate the wearer's own voice from simultaneous ambient sounds. Wearers can speak comfortably and confidently without being bothered by the sound of their own voice or distracted by ambient sounds.
- HandsFree Mode:** By using Bluetooth®, the wearer can use their hearing aids for hands-free calls with their iPhone in different environments. In noisy environments, background noise is reduced, allowing the hearing aid wearer to hear clearly what is being said on the phone. The hearing aid's built-in microphone lets conversation partners hear the wearer's voice in high quality.
- Available in our **Level 7 Signia AX** hearing aids, **Auto EchoShield** creates a cleaner, more immersive environment by neutralizing reverberation and sudden loud noises instantly, before they become an issue in atriums, hallways, lecture halls, and other complex environments.

*Visit signia-pro.com/en-us/claims-data for more details

Performance Levels

7/5/3
and sDemo



24/7 Signia Assistant



Bluetooth Connectivity



Li-ion Rechargeability



Own Voice Processing 2.0



Dynamic Soundscape Processing 2.0



Acoustic Motion Sensors



Signia App



Tinnitus Relief



T-coil



Auto EchoShield



HandsFree for iOS

Available Colors:



Signia app is equipped with the Signia Assistant for 24/7 support plus Signia Face

Mask Mode Improves audibility for those communicating with someone wearing a face mask – all with the touch of a button in the Signia app.

Type	Pure Charge&Go T AX
Available Performance Level	7/5/3
Matrix	S (110/46) M (119/60) P (122/65) HP (131/75)
Battery Size	Lithium-ion power cell
Features	
Frequency channels/ compression controls	7 AX 96/20 (48 channels per processor) 5 AX 64/16 (32 channels per processor) 3 AX 48/12 (24 channels per processor)
Hearing Programs	6
Bluetooth® Streaming	•
TwinPhone ¹	•
TeleCare	•
Tinnitus therapies	•
CROS/BiCROS ²	•
Rocker Switch ³	•
Rechargeable	•
T-coil	•
OVP 2.0	•
Eco-Shield	•
Hands Free	•
App/Accessory Compatibility	
Signia App	•
StreamLine TV/StreamLine Mic	•
miniPocket	•
Noahlink Wireless	•
Fitting set	EarWear 3.0

Pure Charge&Go T AX works exclusively with the new EarWear 3.0.

Multiple charger options available



Pure Charger
Standard Option for Level
3 AX Hearing Aids



Pure Portable Charger
Upgrade Option for Level
5 AX and 7 AX Hearing Aids



Pure Dry&Clean Charger
Upgrade Option for Level 3 AX,
5 AX and 7 AX Hearing Aids

Compatible with:



StreamLine™ TV



StreamLine™ Mic



miniPocket™

¹ Activate via rocker switch, app or accessory

² Compatible with CROS Pure Charge&Go T AX

³ Programmable via Connexx

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners.

- Standard
- Optional



Apple App Store is a servicemark of Apple Inc.

iPhone is a trademark of Apple Inc., registered in the U.S. and other countries.

Android and Google Play are trademarks of Google Inc.

Signia is a registered trademark of WSAUD A/S. 1.0 SIV-13035-22

**Be
Brilliant™**