

# Signia CROS Solutions:

World's-first innovations available to those with single-sided deafness.



signia

## Rechargeable CROS/BiCROS solution with portable charger

- Compatible with:  
Pure Charge&Go AX,  
Pure Charge&Go T AX



## RECHARGEABLE — CROS Pure Charge&Go AX Transmitter

- New chip with increased processing speed for improved synchronization of wireless data exchange
- Own Voice Processing (OVP) in BiCROS fittings
- Directional pattern skewed toward the unaidable ear preserving CROS benefit even in noise
- CROS mixing-level adjustment available in the Signia App to balance spatial benefit vs noise level
- Galvanic rechargeability with on-the-go charging option
- Motion sensors provide improved speech understanding even while moving

## The most flexible rechargeable CROS/BiCROS experience

- Li-ion rechargeability and auto on/off
- Compatible with:  
Pure Charge&Go X,  
Motion Charge&Go S/P/SP X  
(will require separate charger)



## RECHARGEABLE — CROS Pure Charge&Go X Transmitter

- Acoustic-motion sensors for the clearest speech understanding that adapts to the wearer's movements in BiCROS fittings
- Own Voice Processing (OVP™) in BiCROS fittings
- Li-ion rechargeability and auto on/auto off function when removed from charger for easy handling in BiCROS fittings
- Directional pattern skewed toward the unaidable ear preserving CROS benefit even in noise
- CROS mixing-level adjustment available in the Signia App to balance spatial benefit vs. noise level
- On-the-go charging options available

## The easy-to-use, direct streaming CROS/BiCROS experience

- Under the Pure 312: For wearers that prefer a battery driven option
- Compatible with: Pure 312 X,  
Motion Charge&Go S/P/SP  
(will require separate charger)



## SIMPLE & CONNECTED — CROS Pure 312 X Transmitter

- Acoustic motion sensors for the clearest speech understanding that adapts to the wearer's movements in BiCROS fittings
- Own Voice Processing (OVP) in BiCROS fittings
- CROS mixing level adjustment available with the Signia App to balance spatial benefit vs. noise level
- Directional pattern skewed toward the unaidable ear for improved speech understanding in noise

## The discreet, ready-to-wear CROS/BiCROS experience

- Compatible with Silk X




## DISCREET & EASY — CROS Silk™ X Transmitter\*

- Natural hearing experience and exceptional sound quality even while using headphones
- Ultra-discreet housing with click-sleeve options for individualized comfort
- Compatible with the Signia App and miniPocket

## ACOUSTIC MOTION SENSORS | DYNAMIC SOUNDSCAPE PROCESSING | OWN VOICE PROCESSING

### LEGEND



Bluetooth connectivity



Rechargeable



Signia app

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice. The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners. Android is a trademark of Google Inc. Copyright © 2021 Signia GmbH. All rights reserved. Signia is a registered trademark of WSAUD A/S. SHS202165 9/21. SIV-11368-21

\* Motion sensor not available in this product.