

Type	iScan II	
Product details	Ear impression scanner	
Dimensions	Approx. L 341x W 320 x H 125 mm	
Weight	Approx. 7 kg	
Technology	Stripe-light scanner with durable high-performance LED	
Number of impressions	1	
Supply voltage	100 - 240 VAC; 50/60 Hz	
Power consumption	Max. 30 W	
Protection category	IP11	
Accuracy	< 20 µm	
Duration of scan	Approx. 60 - 70 sec. when following recommended computer performance	
Output data format	ASCII or STL	
Ports	USB 2.0; USB 3.0	
Scope of delivery	Scanner, power cable, USB cable, CD with iScan II software, calibration data and user manual, 1 impression holder, printed version of Quick Installation Guide	
System requirements	Minimum	Recommended
System	Windows 7	Windows 10
Ram	2 GB RAM	4 GB RAM
USB	USB Port 2.0	USB Port 3.0
Graphic card	Graphic card minimum 256 MB RAM	Graphic card 512 MB RAM